Instruction Manual Incubating Orbital Shakers

EN - English		1
FR - Français		15
ES - Español		31
IT - Italiano		47
DE - Deutsch		63
PT - Português		79
NL - Nederlands		95
NO - Norsk		99
DA - Dansk		103
SV - Svenska		107
FI - Suomi		111
HU - Magyar		115
PL - Polski		119

European Catalogue Numbers:

Large Capacity Shaker

EU/UK/CH Plug: 444-7084

Model 5000I Incubating Orbital Shaker

EU/UK/CH Plug: 444-7109

Model 5000IR Incubating/Refrigerating Orbital Shaker

EU/UK/CH Plug: 444-7112

TABLE OF CONTENTS

Package Contents	1
Warranty	1
Installation	2
Maintenance & Servicing	2
Intended Use	2
Equipment Disposal	2
Environmental Conditions	2
Safety Instructions	3
Warning	3
Standards & Regulations	3
Control Panel	4
Specifications	5-7
Operating Instructions	8-11
Technical Service	11
Troubleshooting	12
Accessories	13-14
EC Declaration of Conformity	123-124

PACKAGE CONTENTS

Incubating Orbital Shaker
Non-skid rubber mat
234cm detachable power cord
Instruction manual

WARRANTY

VWR International warrants that this product will be free from defects in material and workmanship for a period of two (2) years from date of delivery. If a defect is present, VWR will, at its option and cost, repair, replace, or refund the purchase price of this product to the customer, provided it is returned during the warranty period. This warranty does not apply if the product has been damaged by accident, abuse, misuse, or misapplication, or from ordinary wear and tear. If the required maintenance and inspection services are not performed according to the manuals and any local regulations, such warranty turns invalid, except to the extent, the defect of the product is not due to such non-performance.

Items being returned must be insured by the customer against possible damage or loss. This warranty shall be limited to the aforementioned remedies. IT IS EXPRESSLY AGREED THAT THIS WARRANTY WILL BE IN LIEU OF ALL WARRANTIES OF FITNESS AND IN LIEU OF THE WARRANTY OF MERCHANTABILITY.

Compliance with local laws and regulations

The customer is responsible for applying for and obtaining the necessary regulatory approvals or other authorizations necessary to run or use the Product in its local environment. VWR will not be held liable for any related omission or for not obtaining the required approval or authorization, unless any refusal is due to a defect of the product.

INSTALLATION

Upon receiving the VWR Incubating Orbital Shaker, check to ensure that no damage has occurred during shipment. It is important that any damage that occurred in transport is detected at the time of unpacking. If you do find such damage the carrier must be notified immediately.

After unpacking, it requires two (2) people to lift the Incubating Orbital Shaker, from the bottom, to place on a level bench or table, away from explosive vapors. Ensure that the surface on which the unit is placed will withstand typical heat produced by the unit. Always place the unit on a sturdy work surface.

The Incubating Orbital Shaker is supplied with a power cord that is inserted into the IEC connector on the back of the unit first, then it can be plugged into a properly grounded outlet. The 230V unit plugs into a 230 volt, 50/60 Hz source.

MAINTENANCE & SERVICING

The Incubating Orbital Shaker is built for long, trouble-free, dependable service. No lubrication or other technical user maintenance is required. However at least every three (3) months you should:

- Unplug the unit.
- Remove any accumulated dirt from the base and tray.
- Check all accessible items to make sure they are properly tightened.

The unit should be given the care normally required for any electrical appliance. Avoid wetting or unnecessary exposure to fumes. Spills should be removed promptly. **DO NOT** use a cleaning agent or solvent on the front panel or lid which is abrasive or harmful to plastics, nor one which is flammable. Always ensure the power is disconnected from the unit prior to any cleaning. If the unit ever requires service, contact your VWR representative.

INTENDED USE

The Incubating Orbital Shaker is designed for a variety of shaking applications.

EQUIPMENT DISPOSAL

This equipment must not be disposed of with unsorted waste. It is your responsibility to correctly dispose of the equipment at life-cycle-end by handing it over to an authorized facility for separate collection and recycling. It is also your responsibility to decontaminate the equipment in case of biological, chemical and/or radiological contamination, so as to protect the persons involved in the disposal and recycling of the equipment from health hazards.

For more information about where you can drop off your waste of equipment, please contact your local dealer from whom you originally purchased this equipment. By doing so, you will help to conserve natural and environmental resources and you will ensure that your equipment is recycled in a manner that protects human health.

ENVIRONMENTAL CONDITIONS

Operating Conditions: Indoor use only.

Temperature:	15 to 32°C
Humidity:	maximum 80% relative humidity, non-condensing
Altitude:	0 to 2000 M above sea level

Non-Operating Storage:

Temperature:	-20 to 65°C
Humidity:	maximum 80% relative humidity, non-condensing

Installation Category II and Pollution Degree 2 in accordance with IEC 664.

SAFETY INSTRUCTIONS

Please read the entire instruction manual before operating the Incubating Orbital Shaker.

WARNING! DO NOT use the Incubating Orbital Shaker in a hazardous atmosphere or with hazardous materials for which the unit was not designed. Also, the user should be aware that the protection provided by the equipment may be impaired if used with accessories not provided or recommended by the manufacturer, or used in a manner not specified by the manufacturer.

Exercise extreme caution when moving the unit, there is a risk of pinching or crushing fingers.

Always operate unit on a level surface for best performance and maximum safety.

DO NOT lift unit by the tray, front bezel or lid.

CAUTION! To avoid electrical shock, completely cut off power to the unit by disconnecting the power cord from the unit or unplug from the wall outlet. Disconnect unit from the power supply prior to maintenance and servicing.

Spills should be removed promptly. **DO NOT** immerse the unit for cleaning.

DO NOT operate the unit if it shows signs of electrical or mechanical damage.

CAUTION! The caution hot indicator light warns that the temperature of the air in the chamber is above 40°C. The light will illuminate and remain lit when the temperature of the air in the chamber reaches approximately 40°C. When the heat is turned off, the caution hot indicator light will stay lit until the temperature of the air in the chamber is less than 40°C.

Earth Ground - Protective Conductor Terminal

Alternating Current

STANDARDS & REGULATIONS

VWR international hereby declares under its sole responsibility that the products comply with the requirements of the following directives and associated standards:

Associated EU directives:


EMC directive	2004/108/EC
LVD directive	2006/95/EC
ROHS directive	2011/65/EU
WEEE directive	2002/96/EC

Safety standards:

EN 61010
EN 61010-2-051
EN 61010-2-010
UL 61010-1
CAN/CSA C22.2 No. 61010-1-04

EMC standards:

EN 45501	EN 61000-4-3
EN 55022A	EN 61000-4-4
EN 6100-3-2	EN 61000-4-5
EN 6100-3-3	EN 61000-4-6
EN 61000-4-1	EN 61000-4-11
EN 61000-4-2	

CONTROL PANEL

The front panel of the Incubating Orbital Shaker contains all the controls and displays needed to operate the unit.

A. Standby button/standby indicator light: The standby indicator light will illuminate when the unit is plugged in. The unit will be in standby mode. Press the standby button to activate the temperature, speed and time functions. The standby indicator light will shut off and the temperature, speed and time displays will illuminate. Press the stand-by button again and the unit will once again be in standby mode.

B. Temperature display: Displays the actual/set-point temperatures in conjunction with the actual/set-point indicator lights. **C.** Up/down arrows for set-point control. On/off button starts/stops the heating function. **D.** The heat indicator light will be illuminated when the unit is heating.

E. Actual indicator light: Illuminates when the temperature displayed is the actual temperature of the air in the chamber.

F. Set-point indicator light: Illuminates when the set-point temperature is displayed.

G. Caution hot indicator light: Illuminates when the air temperature of the chamber is above 40°C (104°F).

H. Speed display: Displays the speed of the shaker. **I.** Up/down arrows for set-point control. On/off button starts/stops shaking function. **J.** The speed indicator light will be illuminated when the unit is shaking.

K. Time display: Displays accumulated time (continuous mode) or how much time is remaining (timed mode). The display range is from 0 to 9999 minutes in one (1) second increments. The display will indicate minutes and seconds until the timer reaches 99 minutes and 59 seconds (99:59), then the display will automatically display minutes up to 9999. **L.** Up/down arrows for set-point control. On/off button starts/stops the time function.

SPECIFICATIONS - LARGE CAPACITY INCUBATING SHAKER

Overall dimensions (L x W x H):	64,8 x 35,6 x 40,6cm
Interior dimensions (L x W x H):	34 x 30,5 x 24,1cm
Tray dimensions (L x W):	33 x 27,9cm
Electrical (50/60 Hz):	230 volts, 5 amps, 450 watts
Fuses:	5mm x 20mm, 5 amp quick acting
Temperature range:	ambient +5°C to 65°C
Temperature uniformity:	±0,5°C at 37°C
Speed range:	15 to 500rpm
Speed accuracy:	
Above 100rpm	±1% of set speed
Below 100rpm	±1rpm
Timer:	1 second to 9999 minutes (increased in 1 second increments)
Orbit:	19mm
Maximum weight capacity:	15,9kg
Controls:	see page 4
Tray material:	aluminum
Ship weight:	37,7kg
European Catalogue Number:	EU/UK/CH plug: 444-7084

SPECIFICATIONS - MODEL 5000I INCUBATING SHAKER

Overall dimensions (L x W x H):	81,5 x 67,6 x 59,7cm
Interior dimensions (L x W x H):	52,3 x 63 x 43,2cm
Tray dimensions (L x W):	45,7 x 45,7cm
Electrical (50/60 Hz):	230 volts, 8 amps, 800 watts
Fuses:	5mm x 20mm, 8 amp quick acting
Temperature range:	ambient +5°C to 65°C
Temperature uniformity:	±0,5°C at 37°C
Speed range:	15 to 500rpm
Speed accuracy:	
Above 100rpm	±1% of set speed
Below 100rpm	±1rpm
Timer:	1 second to 9999 minutes (increased in 1 second increments)
Orbit:	25mm
Maximum weight capacity:	22,7kg
Controls:	see page 4
Tray material:	aluminum
Ship weight:	97,5kg
European Catalogue Number:	EU/UK/CH plug: 444-7109

SPECIFICATIONS - MODEL 5000IR INCUBATING/REFRIGERATING SHAKER

Overall dimensions (L x W x H):	104,4 X 67,6 X 59,7cm
Interior dimensions (L x W x H):	52,3 x 63 x 43,2cm
Tray dimensions (L x W):	45,7 x 45,7cm
Electrical (50/60 Hz):	230 volts, 10 amps, 800 watts
Fuses:	5mm x 20mm, 10 amp quick acting
Temperature range:	ambient -15°C to 65°C
Temperature uniformity:	±0,5°C at 37°C
Speed range:	15 to 500rpm
Speed accuracy:	
Above 100rpm	±1% of set speed
Below 100rpm	±1rpm
Timer:	1 second to 9999 minutes (increased in 1 second increments)
Orbit:	25mm
Maximum weight capacity:	22,7kg
Controls:	see page 4
Tray material:	aluminum
Ship weight:	117,9kg
European Catalogue Number:	EU/UK/CH plug: 444-7112

5000IR UPPER AMBIENT CONTROL LIMIT

Please follow the procedure below if setting the 5000IR unit's temperature at/near the room's ambient temperature.

For optimum temperature control in the chamber of the 5000IR unit, the upper ambient control limit is to be set a minimum of 3°C above the measured room temperature.

Example: Measured (ambient) temperature is 22°C, the upper ambient control limit should be set at a minimum of 25°C.

With the unit at the heat/cool setting as described below, and in cool mode ("cool" on time display), a temperature setting is displayed in the temperature window. The default temperature setting of the 5000IR is 28°C. This is the factory set upper ambient control limit. This setting can be adjusted between 20°C and 32°C by the end user.

- Any unit set point above the upper limit setting will heat only
- Any unit set point below 20°C will cool only
- Any set point between 20°C and the selected upper limit setting, up to 32°C, the chiller and heater will both control chamber temperature

Heat/Cool Setting Procedure to Adjust Upper Ambient Control Limit

1. Put unit in standby mode.
2. Press and hold the speed down and time down time arrows simultaneously, now press the standby button. Release all buttons.
3. Time display will not read "cool" and the temperature display will read 28°C. If the time display reads "heat" press the time up/down arrow until the word "cool" appears in the display.
4. Now you can set the upper ambient control limit. Using the temperature up/down arrows adjust this temperature to be 3°C above the measured ambient temperature.
5. Once the desired temperature is selected, press the standby button to return to normal operation.
6. This procedure does not require recalibrating the unit.

Additional 5000IR Notes:

- 5000IR units have a toggle switch on the left side. This switch must be in the "on" position (depress "I") for the unit to function properly.
- For 5000IR units, make sure the chiller's hose is positioned to drain freely (into a sink when possible).
- Do not allow the end of the chiller hose to become submerged.
- If chiller hose is positioned in a drain container, the container should have a volume of at least 2 gallons if the unit will be unattended for 48 hours.

OPERATING INSTRUCTIONS

The Incubating Orbital Shaker has been designed for the temperature, speed and time functions to work independently of one another. The temperature and speed can be reset without resetting the timer and the timer can be stopped and started without interrupting the heating and shaking functions.

1. Getting ready:

- Plug the power cord into a properly grounded outlet. The standby indicator light will illuminate, verifying power to the unit.
- Press the standby button to move the unit from standby mode. The standby indicator light will turn off and the temperature, speed and time displays will illuminate, displaying the previously used settings.

Note: 5000IR units have a toggle switch on the left side. This switch must be in the "on" position (depress "I") for the unit to function properly.

Note: See the procedure on page 8 for 5000IR units that will have the chamber temperature set at/near the ambient temperature.

2. Setting temperature:

- Press the up/down arrows below the temperature display until you reach the desired temperature. When you release the button, the display will blink off and then on indicating the new set temperature has been accepted.
- Press the on/off button to start the heating function. The indicator light below the temperature display will illuminate to indicate the heating function is in use and remain lit until heating has ceased.
- Temperature adjustments can be made without interrupting heating by using the up/down arrows below the temperature display. After the change has been made and you release the button, the display will blink off and then on indicating the new set temperature has been accepted.
- To stop the heating function, press the on/off button below the temperature display. The heat indicator light will turn off.

CAUTION HOT indicator:

The caution hot indicator light warns that the temperature of the air in the chamber is above 40°C. The light will illuminate and remain lit when the temperature of the air in the chamber reaches approximately 40°C. When the heat is turned off, the caution hot indicator light will stay lit until the temperature of the air in the chamber is less than 40°C.

3. Setting speed:

- Press the up/down arrows below the speed display until you reach the desired speed. When you release the button, the display will blink off and then on indicating the new set speed has been accepted.
- Press the on/off button to start the shaking function. The indicator light below the speed display will illuminate and blink until the set-point is reached. Once the set-point is reached the light will stop blinking and remain lit until shaking has ceased. The micro-processor controlled ramping feature slowly increases speed until the set-point is reached which helps to avoid splashing, and provides excellent low end control.
- Speed adjustments can be made without interrupting shaking by using the up/down arrows below the speed display. After the change has been made and you release the button, the display will blink off and then on indicating the new set speed has been accepted.
- To stop the shaking function, press the on/off button below the speed display. The speed indicator light will turn off.

OPERATING INSTRUCTIONS (CONT'D)

4. Setting time to zero (0:00) and continuous mode: Accumulated time.

- a. Press and hold the on/off button below the time display. After three (3) seconds the display will indicate the previous set time.
- b. Simultaneously press both the up and the down arrows, the display will indicate zero (0:00). The unit time is now set to zero (0:00) minutes. Alternately, you can use the up/down arrows to get to zero (0:00).
- c. Press the on/off button below the time display. The display will indicate accumulated time. The up/down arrows will become inactive. To stop timer, press the on/off button again. **IMPORTANT:** This will **NOT** interrupt the shaking function. Press the on/off button below the speed display to interrupt the shaking function.
- d. To reset, press and hold the on/off button below the time display. After three (3) seconds the display will indicate the previous set time, which was zero (0:00).

5. Setting timed mode: Programmed time.

- a. Press the up/down arrows below the time display until you reach the desired time.
- b. Start this function by pressing the on/off button below the time display. The unit will run for the selected time, the up/down arrows will become inactive while the timer is running. The unit will stop shaking when time display reaches zero (0:00). Four (4) audible beeps will indicate the count down function is complete. The time display will default back to the set time. To repeat for the same time, simply press the on/off button again.
- c. To interrupt an automatic timing cycle before it is completed, press the on/off button below the time display. The display will flash off and on to indicate the time function is on "hold". **IMPORTANT:** This will **NOT** interrupt the shaking function. Press the on/off button below the speed display to interrupt the shaking function. Restart the timer by pressing the on/off button below the time display. Unit will continue counting down to zero (0:00). When the display reaches zero (0:00), you will hear the four (4) audible beeps that indicate the count down function is complete and the shaking function will cease.

6. Turning unit off:

- a. To turn the unit off, press the standby button. The temperature, speed and time displays will be blank, the standby indicator light will illuminate. The Incubating Orbital Shaker should be kept in standby mode when not in use. To completely cut off power to the unit, disconnect the power cord from the unit or unplug from the wall outlet.

OPERATING TIPS

Opening the lid on the Incubating Orbital Shaker will cause the unit to pause shaking and/or heating. Close the lid and the unit will resume shaking and/or heating at current settings.

Centering your sample and even weight distribution on the tray helps with balance and stability.

As a safety feature, a built-in program will shut power off to the motor if the tray is prevented from rotating, or the unit is overloaded beyond its recommended weight capacity.

The shaker will automatically restart after a power interruption. Built-in memory maintains the last used temperature, speed and time settings during a power interruption.

LOAD SENSING FUNCTION

The Incubating Orbital Shaker is equipped with a load sensing function that can be activated by the user. This function provides protection against improper positioning of load and maximum load being exceeded. When activated, the unit will automatically sense improper load conditions and slow to a safe running speed, then display that speed followed by an E04 error message on the speed display. The unit will also beep three (3) times every 60 seconds until the error is reset by pressing the speed on/off button. To activate the load sensing function use the following steps:

1. Place the unit in standby mode.
2. Press and hold the speed on/off button and press the standby button. The unit will beep two (2) times, confirming the function is enabled.

OPERATING INSTRUCTIONS (CONT'D)

3. To restore normal operation, remove AC power to the unit for ten (10) seconds and then restore.

If the E04 error occurs be sure the load is within the maximum specification and properly balanced (centered on tray) and/or reduce sample size/speed before restarting the unit. If the E04 error occurs due to acceptable sample vibration or another vibration source, the vibration sensing function can be disabled as described above.

BEEPER PREFERENCE (muting audible alarm)

To silence beeper operation (except for error codes), with the unit in standby mode, press and hold the time on/off button and press the standby button. Release the standby button first, and then release the on/off button. To restore normal beeper operation, remove AC power to unit for ten (10) seconds and then restore.

SINGLE POINT CALIBRATION PROCEDURE

This procedure is used to fine tune and calibrate the Incubating Orbital Shaker at a specific temperature setting. This process may be repeated for up to three (3) separate set-points. If a fourth calibration set-point is entered, the first set-point entered will be overwritten.

1. Turn unit on.
2. Set desired temperature.
3. Stabilize one (1) hour or more, measuring the temperature with a temperature probe/thermometer.
4. Press and hold standby button, then press the temperature up arrow once. The unit will beep two (2) times, confirming calibration mode. The display will now be flashing.
5. Press the temperature up/down arrows until the display matches the temperature probe/thermometer.
6. Press standby button to exit calibration mode and return to normal heating.

This process may be repeated at the same set-point, multiple times for fine tuning if desired.

The unit will now use the biased offset for that specific temperature setting and increase or decrease temperature accordingly to bring the temperature to set temperature. The decimal point of the display will flash to indicate a biased offset is being used. All other temperature settings will use the standard internal calibration. This offset will be stored in memory and retained until reset.

To restore unit to factory setting:

Press and hold the standby button while pressing the temperature down arrow once. The reset will be confirmed with two (2) audible beeps. Press the standby button to exit calibration mode and return to normal heating.

Speed Calibration Procedure:

This procedure is used to self calibrate the incubating shakers. The tray should be free of any samples, vessels, and accessories prior to calibrating.

1. Turn unit on. Speed and time displays will be illuminated.
2. Press and hold the standby button and momentarily press the speed on/off button. The speed display should read "CAL".
3. The unit will be run for approximately one (1) minute and automatically calibrate.

RS-232 SERIAL PORT

The RS-232 serial port provides two-way communications for data logging and unit control.

TECHNICAL SERVICE

For more information or technical assistance contact your local VWR representative or visit www.vwr.com.

Web resources: Visit VWR's website at www.vwr.com for:

- Complete technical service contact information
- Access to VWR's Online Catalogue, and information about accessories and related products
- Additional product information and special offers

TROUBLESHOOTING

To clear error press the standby button


Problem	Cause	Solution
Unit fails to power on	Missing or blown fuse	Check and replace fuse if necessary If problem persists, please contact your VWR representative for repair.
Unit is excessively noisy	Rattling or ticking sounds may indicate a loose screw on the tray	Ensure that tray is secured tightly. If problem persists, please contact your VWR representative for repair.
Unit not shaking at proper speed	-	Perform speed calibration on page 11 If problem persists, please contact your VWR representative for repair.
Unit not heating to proper temperature	-	Perform single point calibration on page 11 If problem persists, please contact your VWR representative for repair.
E1	RTD open or $\pm 1^{\circ}\text{C}$ temperature deviation from set-point (after unit has stabilized)	This error cannot be fixed by the end user. If problem persists, please contact your VWR representative for repair.
E2	RTD shorted or temperature below 0°C	This error cannot be fixed by the end user. If problem persists, please contact your VWR representative for repair.
E3	Mechanical obstruction Drive system failure Ceased bearing Drive belt broken	Remove mechanical obstruction. If problem persists, the reason may be the drive system. Please contact your VWR representative for repair.
E4	Improper positioning of load or Maximum load exceeded	Ensure the load is evenly distributed and does not exceed the maximum load capacity for the unit. See "Load Sensing Function" on page 10-11. Verify load has not exceeded the maximum allowable weight. Reduce load if necessary. If error code resumes, contact your VWR representative for repair.
E6	Over temperature error	This error cannot be fixed by the end user. Please contact your VWR representative for repair.
E7	Default sensor Error (5000IR units only)	Error will shut unit off if set temperature is 20°C or below. Unit is still operational with a set temperature above 20°C . Error will be indicated every 5 minutes until repaired. Please contact your VWR representative for repair.

ACCESSORIES

PLATFORMS

Platform can be used for mounting Erlenmeyer flask/media bottle clamps and test tube racks. Slides over-top included tray. Tighten with adjustment screws. Stainless steel construction.

DESCRIPTION	DIMENSIONS (L x W)	CAT. No.
Platform	33 x 28cm	444-7087
Platform	47,5 x 47,5cm	444-0697

TEST TUBE RACKS - Requires Platform accessory.

		33 x 28cm	47,5 x 47,5	
DESCRIPTION	TUBE CAPACITY	PLATFORM CAPACITY	PLATFORM CAPACITY	CAT. No.
Half Size Test Tube Racks				
10 to 13mm Test Tube Rack	63	2	4	444-0278
14 to 16mm Test Tube Rack	48	2	4	444-0279
18 to 20mm Test Tube Rack	35	2	4	444-0280
22 to 25mm Test Tube Rack	24	2	4	444-0281
15mL Centrifuge Tube Rack	35	2	4	444-0283
50mL Centrifuge Tube Rack	12	2	4	444-0284
1,5 to 2mL Micro-Tube Rack	70	2	4	444-0276
Full Size Stationary Test Tube Racks				
10 to 14mm Test Tube Rack	48	N/A	3	444-7037
16 to 20mm Test Tube Rack	33	N/A	3	444-7038
21 to 25mm Test Tube Rack	21	N/A	3	444-7039
Full Size Pivoting Test Tube Racks				
13mm Test Tube Rack	90	N/A	2	444-2955
16mm Test Tube Rack	60	N/A	2	444-2956
20mm Test Tube Rack	40	N/A	2	444-2957
25mm Test Tube Rack	24	N/A	2	444-2958
30mm Test Tube Rack	21	N/A	2	444-2959

Installation of Half Size Test Tube Racks:

1. Unit should be in standby mode.
2. Remove non-skid rubber mat. Slide Platform over-top included tray.
3. Using the hardware provided, attach brackets to the Platform. Position brackets as shown. The proper spacing between the brackets is 19,7cm.
4. Align the test tube rack side-to-side with the brackets on the Platform.
5. Press the test tube rack in place, you will hear the test tube rack clip into position. No mounting screws are required.

Installation of Full Size Pivoting Test Tube Racks (Model 5000I/5000IR Only):

1. Unit should be in standby mode.
2. Line up the 4 openings on the bottom of the outer rack with the four compatible openings on the Universal Platform. Secure with the four screws provided.
3. If a tilt angle is desired, loosen wing knobs, tilt inner holder and tighten knobs.
4. To maximize the platforms capacity for holding racks in the stationary position, remove the outer holder and secure the inner holder to the Universal Platform using the four openings on the bottom of the inner holder.

Installation of Full Size Stationary Test Tube Racks (Model 5000I/5000IR Only):


1. Line up the two openings on the sides of the test tube rack with the two compatible openings on the Universal Platform.
2. Place the washer on top of the opening on the test tube rack, place the clamp knob through the washer and rack and hand tighten the clamp knob to the platform.

ACCESSORIES

PLATFORMS

Platform can be used for mounting Erlenmeyer flask/media bottle clamps and test tube racks. Slides over-top included tray. Tighten with adjustment screws. Stainless steel construction.

DESCRIPTION	DIMENSIONS (L x W)	CAT. No.
Platform	33 x 28cm	444-7087
Platform	47,5 x 47,5cm	444-0697

FLASK CLAMPS - Requires Platform accessory

33 X 28cm		47,5 x 47,5		33 X 28cm		47,5 x 47,5	
DESCRIPTION	PLATFORM CAPACITY	PLATFORM CAPACITY	CAT. No.	DESCRIPTION	PLATFORM CAPACITY	PLATFORM CAPACITY	CAT. No.
Stainless Steel Flask Clamps				Stainless Steel Media Bottle Clamps			
10mL Erlenmeyer Flask Clamp	60	113	444-7040	500mL Media Bottle Clamp	5	16	444-7035
25mL Erlenmeyer Flask Clamp	25	64	444-7041	1L Media Bottle Clamp	2	10	444-7036
50mL Erlenmeyer Flask Clamp	13	32	444-7042	PVC Flask Clamps			
125mL Erlenmeyer Flask Clamp	10	20	444-7043	125mL Erlenmeyer Flask Clamp	10	20	444-2960
250mL Erlenmeyer Flask Clamp	9	20	444-7044	250mL Erlenmeyer Flask Clamp	8	18	444-2961
500mL Erlenmeyer Flask Clamp	7	13	444-7045	500mL Erlenmeyer Flask Clamp	5	12	444-2962
1L Erlenmeyer Flask Clamp	4	8	444-7030	1L Erlenmeyer Flask Clamp	2	8	444-2963
2L Erlenmeyer Flask Clamp	N/A	5	444-7031	2L Erlenmeyer Flask Clamp	N/A	4	444-2964
2.8L Erlenmeyer Flask Clamp	N/A	2	444-7032				
4L Erlenmeyer Flask Clamp	N/A	4	444-7033				
6L Erlenmeyer Flask Clamp	N/A	2	444-7034				

Installation of Flask Clamps:

1. Unit should be in standby mode.
2. Remove non-skid rubber mat. Slide Platform over-top included tray.
3. Line up the hole(s) in the flask clamp to the hole(s) in the Platform. Using the screw(s) provided, hand tighten into place using a flathead screwdriver. **DO NOT** over tighten.
4. Insert flask/media bottle into clamp and wrap the spring around the neck of the clamp to hold secure.

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hereby declare that the products, mentioned in Schedule comply with the essential requirements of the following Directives of the European Parliament and council:

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Title of Council Directive on harmonization of the laws of the member states	
Directive 2006/95/EC of the European Parliament and of the Council of 12 December 2006 on the harmonisation of the laws of Member States relating to Electrical Equipment designed for use within certain voltage limits	
Directive 2004/108/EC of the European Parliament and of the Council of 15 December 2004 on the approximation of the laws of the Member States relating to electromagnetic compatibility	
and are in conformity with the following standard(s) and/or other normative document(s): <i>et sont conformes à la (aux) norme(s) ou autre(s) document(s) normatif(s) :</i> (den) folgenden Normen oder anderer normativer Dokumente entsprechen:	
Reference number	Title of the applied standards:
IEC 61010-2-010	Safety requirements for electrical equipment for measurement, control and laboratory use – Part 2-010: Particular requirements for laboratory equipment for heating of materials
IEC 61010-1	Safety requirements for electrical equipment for measurement, control and laboratory use – Part 1: General requirements – Second Edition
IEC 61010-2-051	Safety requirements for electrical equipment for measurement, control and laboratory use – Part 2-051: Particular requirements for laboratory equipment for mixing and stirring
IEC 61000-3-2; -3-3 IEC 61000-4-2; -4-3; -4-4 IEC 61000-4-5; -4-6; -4-11 EN 45501; 55022A	Electromagnetic compatibility (EMC)

Year of CE marking: 2006

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Place and date of issue: Leuven, 11/05/2012

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Yves Van Damme
 European Regulatory Affairs Manager

Frankie Vangeel
 Managing Director

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Schedule 1
Annexe 1
Anhang 1

Article Number	Article description
444-7084	Shaker Incubated Large Capacity, Digital, 230V, EU
444-7085	Shaker Incubated Large Capacity, Digital, 230V, UK
444-7086	Shaker Incubated Large Capacity, Digital, 230V, CH

2/2

VWR International bvba/sprl Researchpark Haasrode 2020, Geldenaaksebaan 464, B-3001 Leuven, Belgium
Tel.: 016 385 011 • Fax: 016 385 385 • E-mail: info@be.vwr.com

Revision 01 of May 2012

EC Declaration of Conformity

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We, **VWR International bvba**,
 Nous, **Researchpark Haasrode 2020**,
 Wir, **Geldenaaksbaan 464**,
B-3001 Leuven

hereby declare that the products, mentioned in schedule 1 comply with the essential requirements of the following Directives of the European Parliament and council:

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Yves Van Damme
European Regulatory Affairs Manager

Kris Nijs
Managing Director

VWR International bvba/srl Researchpark Haasrode 2020, Geldenaaksebaan 464, B-3001 Leuven, Belgium
Tel.: 016 385 011 • Fax: 016 385 385 • E-mail: info@be.vwr.com

2/3

Revision 01 of April 2014

EC Declaration of Conformity

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Schedule 1
Annexe 1
Anhang 1

Article Number	Article description
444-7109	SHAKER INCUBATING MODEL 5000i 230V EU SHAKER INCUBATING MODEL 5000i 230V EU SHAKER INCUBATING MODEL 5000i 230V EU
444-7110	SHAKER INCUBATING MODEL 5000i 230V UK SHAKER INCUBATING MODEL 5000i 230V UK SHAKER INCUBATING MODEL 5000i 230V UK
444-7111	SHAKER INCUBATING MODEL 5000i 230V CH SHAKER INCUBATING MODEL 5000i 230V CH SHAKER INCUBATING MODEL 5000i 230V CH

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

Yves Van Damme
European Regulatory Affairs Manager

Kris Nijs
Managing Director

VWR International bvba/sprl Researchpark Haasrode 2020, Geldenaaksebaan 464, B-3001 Leuven, Belgium
 Tel.: 016 385 011 • Fax: 016 385 385 • E-mail: info@bs.vwr.com

Revision 01 of April 2014

EC Declaration of Conformity

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Schedule 1
Annexe 1
Anhang 1

Article Number	Article description
444-7112	SHAKER INC/REFRIG MODEL 5000IR 230V EU SHAKER INC/REFRIG MODEL 5000IR 230V EU SHAKER INC/REFRIG MODEL 5000IR 230V EU
444-7113	SHAKER INC/REFRIG MODEL 5000IR 230V UK SHAKER INC/REFRIG MODEL 5000IR 230V UK SHAKER INC/REFRIG MODEL 5000IR 230V UK
444-7114	SHAKER INC/REFRIG MODEL 5000IR 230V CH SHAKER INC/REFRIG MODEL 5000IR 230V CH SHAKER INC/REFRIG MODEL 5000IR 230V CH

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Ireland / Northern Ireland

VWR International Ltd /
VWR International (Northern Ireland) Ltd
Orion Business Campus
Northwest Business Park
Ballycoolin
Dublin 15
Tel.: 01 88 22 222
Fax: 01 88 22 333
Email: sales@ie.vwr.com

Italy

VWR International PBI S.r.l.
Via San Giusto 85
20153 Milano (MI)
Tel.: 02-3320311/02-487791
Fax: 800 152999/02-40090010
Email: info@it.vwr.com

The Netherlands

VWR International B.V.
Postbus 8198
1005 AD Amsterdam
Tel.: 020 4808 400
Fax: 020 4808 480
Email: info@nl.vwr.com

Norway

VWR International AS
Haavard Martinsens vei 30
0978 Oslo
Tel.: 02290
Fax: 815 00 940
Email: info@no.vwr.com

Poland

VWR International Sp. z o.o.
Limbowa 5
80-175 Gdansk
Tel.: 058 32 38 200 do 204
Fax: 058 32 38 205
Email: info@pl.vwr.com

Portugal

VWR International -
Material de Laboratório, Lda
Edifício Neopark
Av. Tomás Ribeiro, 43- 3 D
2790-221 Carnaxide
Tel.: 21 3600 770
Fax: 21 3600 798/9
Email: info@pt.vwr.com

Spain

VWR International Eurolab S.L.
C/ Tecnología 5-17
A-7 Linars Park
08450 - Llinars del Vallès
Barcelona
Tel.: 902 222 897
Fax: 902 430 657
Email: info@es.vwr.com

Sweden

VWR International AB
Fagerstagatan 18a
163 94 Stockholm
Tel.: 08 621 34 00
Fax: 08 621 34 66
Email: kundservice@se.vwr.com

Switzerland

VWR International GmbH
Lerzenstrasse 16/18
8953 Dietikon
Tel.: 044 745 13 13
Fax: 044 745 13 10
Email: info@ch.vwr.com

Turkey

VWR International Laboratuvar
Teknolojileri Ltd.İti.
Orta Mah. Cemal Gürsel Caddesi
Ördekcioglu İlmerkezi No.32/1
34896 Pendik - Istanbul
Tel.: +90216 598 2900
Fax: +90216 598 2907
Email: info@pro-lab.com.tr

UK

VWR International Ltd
Customer Service Centre
Hunter Boulevard - Magna Park
Lutterworth
Leicestershire
LE17 4XN
Tel.: 0800 22 33 44
Fax: 01455 55 85 86
Email: uksales@uk.vwr.com

Australia

VWR International, Pty Ltd.
Level 1, Unit 1a/60 Enterprise Place
Tingalpa, Queensland, 4173
Tel.: 1300 727 696
Fax: 1300 135 123

China

VWR International China Co., Ltd
Rm.219, 2100 Dongming Road
Pudong New District
Shanghai 200123
Tel.: +86-21-5898 6888
Fax: +86-21-5855 8801
Email: info_china@vwr.com

India

VWR Lab Products Private Limited
No.139, BDA Industrial Suburb,
6th Main, Tumkur Road, Peenya
Post,
Bangalore, India – 560058
Tel.: +91-80-28078400
Fax: +91-80-28078410
Email: vwr_india@vwr.com

New Zealand

Global Science - A VWR Company
241 Bush Road
Albany 0632, Auckland
Tel.: 0800 734 100
Fax: 0800 999 002
Email: sales@globalscience.co.nz

Singapore

VWR Singapore Pte Ltd
18 Gul Drive
Singapore 629468
Tel: +65 6505 0760
Fax: +65 6264 3780
Email: sales@sg.vwr.com

Legal Manufacturer:

VWR International bvba • Researchpark Haasrode 2020 • Geldenaaksebaan 464 • B-3001 Leuven • + 32 16 385011 • <http://be.vwr.com>
715085-00 (Rev 6)