Incubating Orbital Shaker

Revision 9/07

Instruction Manual

Euro plug: 444-7084

UK plug: 444-7085

Swiss plug: 444-7086

CE

TABLE OF CONTENTS

Package Contents	1
Installation	1
Maintenance & Servicing	2
Environmental Conditions	2
Intended Use	2
Equipment Disposal	2
Safety Instructions	3
Standards & Regulations	3
Control Panel	4
Specifications	5
Operating Instructions	5-8
Troubleshooting	8
Replacement Parts	9-10
Accessories	11-12

PACKAGE CONTENTS

Incubating Orbital Shaker
Non-skid rubber mat
234cm detachable power cord
Instruction manual

INSTALLATION

Upon receiving the VWR Incubating Orbital Shaker, check to ensure that no damage has occurred during shipment. It is important that any damage that occurred in transport is detected at the time of unpacking. If you do find such damage the carrier must be notified immediately.

After unpacking, it requires two (2) people to lift the Incubating Orbital Shaker, from the bottom, to place on a level bench or table, away from explosive vapors. Ensure that the surface on which the unit is placed will withstand typical heat produced by the unit. Always place the unit on a sturdy work surface.

The Incubating Orbital Shaker is supplied with a power cord that is inserted into the IEC connector on the back of the unit first, then it can be plugged into a properly grounded outlet. The 230V unit plugs into a 230 volt, 50/60 Hz source.

MAINTENANCE & SERVICING

The Incubating Orbital Shaker is built for long, trouble-free, dependable service. No lubrication or other technical user maintenance is required. However at least every three (3) months you should:

- Unplug the unit.
- Remove any accumulated dirt from the base and tray.
- Check all accessible items to make sure they are properly tightened.

The unit should be given the care normally required for any electrical appliance. Avoid wetting or unnecessary exposure to fumes. Spills should be removed promptly. **DO NOT** use a cleaning agent or solvent on the front panel or lid which is abrasive or harmful to plastics, nor one which is flammable. Always ensure the power is disconnected from the unit prior to any cleaning. If the unit ever requires service, contact your VWR representative.

ENVIRONMENTAL CONDITIONS

Operating Conditions: Indoor use only.

Temperature: 15 to 32°C (59 to 90°F)

Humidity: maximum 80% relative humidity, non-condensing

Altitude: 0 to 6,562 ft (2000 M) above sea level

Non-Operating Storage:

Temperature: -20 to 65°C (-4 to 149°F)

Humidity: maximum 80% relative humidity, non-condensing

Installation Category II and Pollution Degree 2 in accordance with IEC 664.

INTENDED USE

The Incubating Orbital Shaker is designed for a variety of shaking applications.

EQUIPMENT DISPOSAL

This equipment must not be disposed of with unsorted waste. It is your responsibility to correctly dispose of the equipment at life-cycle-end by handing it over to an authorized facility for separate collection and recycling. It is also your responsibility to decontaminate the equipment in case of biological, chemical and/or radiological contamination, so as to protect the persons involved in the disposal and recycling of the equipment from health hazards.

For more information about where you can drop off your waste of equipment, please contact your local dealer from whom you originally purchased this equipment. By doing so, you will help to conserve natural and environmental resources and you will ensure that your equipment is recycled in a manner that protects human health.

SAFETY INSTRUCTIONS

Please read the entire instruction manual before operating the Incubating Orbital Shaker.

WARNING! DO NOT use the Incubating Orbital Shaker in a hazardous atmosphere or with hazardous materials for which the unit was not designed. Also, the user should be aware that the protection provided by the equipment may be impaired if used with accessories not provided or recommended by the manufacturer, or used in a manner not specified by the manufacturer.

Always operate unit on a level surface for best performance and maximum safety.

DO NOT lift unit by the tray or lid.

CAUTION! To avoid electrical shock, completely cut off power to the unit by disconnecting the power cord from the unit or unplug from the wall outlet. Disconnect unit from the power supply prior to maintenance and servicing.

Spills should be removed promptly. **DO NOT** immerse the unit for cleaning.

DO NOT operate the unit if it shows signs of electrical or mechanical damage.

Earth Ground - Protective Conductor Terminal

Alternating Current

STANDARDS & REGULATIONS

VWR International hereby declares under its sole responsibility that the construction of this product conforms in accordance with the following standards:

Safety standards:

IEC 61010-1	Safety requirements for electrical equipment for measurement, control and laboratory use. Part I: General Requirements.
IEC 61010-2-010	Part II: Particular requirements for laboratory equipment for the heating of materials.
IEC 61010-2-051	Part II: Particular requirements for laboratory equipment for mixing and stirring.
UL Std. No. 61010-1	
CSA/CAN C22.2 No. 0-M91	
CSA/CAN C22.2 No. 61010-1-04	


EMC standards:

FCC-B	EN45501
EN55022-A	EN6100-3-3/3-3
EN61000-4-2	EN61000-4-3
EN61000-4-4	EN61000-4-5
EN61000-4-6	EN61000-4-11

Associated EU guidelines:

EMC directive 89/336/EEC
LVD directive 73/23/EEC

Patent Pending

CONTROL PANEL

The front panel of the Incubating Orbital Shaker contains all the controls and displays needed to operate the unit.

A. Standby button/standby indicator light: The standby indicator light will illuminate when the unit is plugged in. The unit will be in standby mode. Press the standby button to activate the temperature, speed and time functions. The standby indicator light will shut off and the temperature, speed and time displays will illuminate. Press the standby button again and the unit will once again be in standby mode.

B. Temperature display: Displays the actual/set-point temperatures in conjunction with the actual/set-point indicator lights. **C.** Up/down arrows for set-point control. On/off button starts/stops the heating function. **D.** The heat indicator light will be illuminated when the unit is heating.

E. Actual indicator light: Illuminates when the temperature displayed is the actual temperature of the air in the chamber.

F. Set-point indicator light: Illuminates when the set-point temperature is displayed.

G. Caution hot indicator light: Illuminates when the air temperature of the chamber is above 40°C (104°F).

H. Speed display: Displays the speed of the shaker. **I.** Up/down arrows for set-point control. On/off button starts/stops shaking function. **J.** The speed indicator light will be illuminated when the unit is shaking.

K. Time display: Displays accumulated time (continuous mode) or how much time is remaining (timed mode). The display range is from 0 to 9,999 minutes in one (1) second increments. The display will indicate minutes and seconds until the timer reaches 99 minutes and 59 seconds (99:59), then the display will automatically display minutes up to 9,999. **L.** Up/down arrows for set-point control. On/off button starts/stops the time function.

Euro plug: 444-7084 UK plug: 444-7085 Swiss plug: 444-7086

SPECIFICATIONS

Overall dimensions (L x W x H):	64.8 x 35.6 x 40.6cm
Interior dimensions (L x W x H):	34 x 30.5 x 24.1cm
Tray dimensions (L x W):	33 x 27.9cm
Electrical (50/60 Hz):	230 volts, 5 amps, 450 watts
Fuses:	5mm x 20mm, 5 amp quick acting
Temperature range:	ambient +5°C to 65°C
Temperature uniformity:	±0.5°C at 37°C
Speed range:	15 to 500rpm
Speed accuracy:	±1%
Timer:	1 second to 9999 minutes (increased in 1 second increments)
Orbit:	19mm
Maximum weight capacity:	15.9kg
Controls:	see page 4
Tray material:	aluminum
Ship weight:	37.7kg

OPERATING INSTRUCTIONS

The Incubating Orbital Shaker has been designed for the temperature, speed and time functions to work independently of one another. The temperature and speed can be reset without resetting the timer and the timer can be stopped and started without interrupting the heating and shaking functions.

1. Getting ready:

- Plug the power cord into a properly grounded outlet. The standby indicator light will illuminate, verifying power to the unit.
- Press the standby button to move the unit from standby mode. The standby indicator light will turn off and the temperature, speed and time displays will illuminate, displaying the previously used settings.

2. Setting temperature:

- Press the up/down arrows below the temperature display until you reach the desired temperature. When you release the button, the display will blink off and then on indicating the new set temperature has been accepted.
- Press the on/off button to start the heating function. The indicator light below the temperature display will illuminate to indicate the heating function is in use and remain lit until heating has ceased.
- Temperature adjustments can be made without interrupting heating by using the up/down arrows below the temperature display. After the change has been made and you release the button, the display will blink off and then on indicating the new set temperature has been accepted.
- To stop the heating function, press the on/off button below the temperature display. The heat indicator light will turn off.

**CAUTION HOT indicator:**

The caution hot indicator light warns that the temperature of the air in the chamber is above 40°C (104°F). The light will illuminate and remain lit when the temperature of the air in the chamber reaches approximately 40°C (104°F). When the heat is turned off, the caution hot indicator light will stay lit until the temperature of the air in the chamber is less than 40°C (104°F).

3. Setting speed:

- a. Press the up/down arrows below the speed display until you reach the desired speed. When you release the button, the display will blink off and then on indicating the new set speed has been accepted.
- b. Press the on/off button to start the shaking function. The indicator light below the speed display will illuminate and blink until the set-point is reached. Once the set-point is reached the light will stop blinking and remain lit until shaking has ceased. The micro-processor controlled ramping feature slowly increases speed until the set-point is reached which helps to avoid splashing, and provides excellent low end control.
- c. Speed adjustments can be made without interrupting shaking by using the up/down arrows below the speed display. After the change has been made and you release the button, the display will blink off and then on indicating the new set speed has been accepted.
- d. To stop the shaking function, press the on/off button below the speed display. The speed indicator light will turn off.

4. Setting time to zero (0:00) and continuous mode: Accumulated time.

- a. Press and hold the on/off button below the time display. After three (3) seconds the display will indicate the previous set time.
- b. Simultaneously press both the up and the down arrows, the display will indicate zero (0:00). The unit time is now set to zero (0:00)

minutes. Alternately, you can use the up/down arrows to get to zero (0:00).

- c. Press the on/off button below the time display. The display will indicate accumulated time. The up/down arrows will become inactive. To stop timer, press the on/off button again. **IMPORTANT:** This will **NOT** interrupt the shaking function. Press the on/off button below the speed display to interrupt the shaking function.
- d. To reset, press and hold the on/off button below the time display. After three (3) seconds the display will indicate the previous set time, which was zero (0:00).

5. Setting timed mode: Programmed time.

- a. Press the up/down arrows below the time display until you reach the desired time.
- b. Start this function by pressing the on/off button below the time display. The unit will run for the selected time, the up/down arrows will become inactive while the timer is running. The unit will stop shaking when time display reaches zero (0:00). Four (4) audible beeps will indicate the count down function is complete. The time display will default back to the set time. To repeat for the same time, simply press the on/off button again.
- c. To interrupt an automatic timing cycle before it is completed, press the on/off button below the time display. The display will flash off and on to indicate the time function is on "hold". **IMPORTANT:** This will **NOT** interrupt the shaking function. Press the on/off button below the speed display to interrupt the shaking function. Restart the timer by pressing the on/off button below the time display. Unit will continue counting down to zero (0:00). When the display reaches zero (0:00), you will hear the four (4) audible beeps that indicate the count down function is complete and the shaking function will cease.

6. Turning unit off:

- a. To turn the unit off, press the standby button. The temperature, speed and time displays will be blank, the standby indicator light will illuminate. The Incubating Orbital Shaker should be kept in standby mode when not in use. To completely cut off power to the unit, disconnect the power cord from the unit or unplug from the wall outlet.

OPERATING TIPS

Opening the lid on the Incubating Orbital Shaker will cause the unit to pause shaking and/or heating. Close the lid and the unit will resume shaking and/or heating at current settings.

Centering your sample and even weight distribution on the tray helps with balance and stability.

As a safety feature, a built-in program will shut power off to the motor if the tray is prevented from rotating, or the unit is overloaded beyond its recommended weight capacity.

The shaker will automatically restart after a power interruption. Built-in memory maintains the last used temperature, speed and time settings during a power interruption.

LOAD SENSING FUNCTION

The Incubating Orbital Shaker is equipped with a load sensing function that can be activated by the user. This function provides protection against improper positioning of load and maximum load being exceeded. When activated, the unit will automatically sense improper load conditions and slow to a safe running speed, then display that speed followed by an E04 error message on the speed display. The unit will also beep three (3) times every 60 seconds until the error is reset by pressing the speed  on/off button. To activate the load sensing function use the

following steps:

1. Place the unit in standby mode.
2. Press and hold the speed on/off button and press the standby button. The unit will beep two (2) times, confirming the function is enabled.
3. To restore normal operation, remove AC power to the unit for ten (10) seconds and then restore.

If the E04 error occurs be sure the load is within the maximum specification and properly balanced (centered on tray) and/or reduce sample size/speed before restarting the unit. If the E04 error occurs due to acceptable sample vibration or another vibration source, the vibration sensing function can be disabled as described above.

BEEPER PREFERENCE (muting audible alarm)

To silence beeper operation (except for error codes), with the unit in standby mode, press and hold the time on/off button and press the standby button. Release the standby button first, and then release the on/off button. To restore normal beeper operation, remove AC power to unit for ten (10) seconds and then restore.

SINGLE POINT CALIBRATION PROCEDURE

This procedure is used to fine tune and calibrate the Incubating Orbital Shaker at a specific temperature setting. This process may be repeated for up to three (3) separate set-points. If a fourth calibration set-point is entered, the first set-point entered will be overwritten.

1. Turn unit on.
2. Set desired temperature.
3. Stabilize one (1) hour or more, measuring the temperature with a temperature probe/thermometer.

4. Press and hold standby button, then press the temperature up arrow once. The unit will beep two (2) times, confirming calibration mode. The display will now be flashing.
5. Press the temperature up/down arrows until the display matches the temperature probe/thermometer.
6. Press standby button to exit calibration mode and return to normal heating.

This process may be repeated at the same set-point, multiple times for fine tuning if desired.

The unit will now use the biased offset for that specific temperature setting and increase or decrease temperature accordingly to bring the temperature to set temperature. The decimal point of the display will flash to indicate a biased offset is being used. All other temperature settings will use the standard internal calibration. This offset will be stored in memory and retained until reset.

To restore unit to factory setting:

Press and hold the standby button while pressing the temperature down arrow once. The reset will be confirmed with two (2) audible beeps. Press the standby button to exit calibration mode and return to normal heating.

RS-232 SERIAL PORT

The RS-232 serial port provides two-way communications for data logging and unit control. If you need additional details, please contact Troemner for technical support at: tech_support@troemner.com.

TROUBLESHOOTING

During operation, any rattling or ticking sounds may indicate a loose screw on the tray, a tray attachment or an accessory. All accessories should be sufficiently tightened in place before starting the unit.

Error Code	Cause
E04	improper positioning of load maximum load exceeded
See 'Load Sensing Function' (page 7). If the E04 error persists, switch the unit off and contact your VWR representative for repairs.	
Error Code	Cause
E03	drive system failure mechanical obstruction

Press the standby button to clear this error and remove the mechanical obstruction. If the E03 error persists the reason may be a drive system failure and should **NOT** be addressed by the end user. Switch the unit off and contact your VWR representative for repairs.

Error Code	Cause
E02	RTD shorted or temperature below 0°C (32°F)
E01	RTD open or ±1°C temperature deviation from set-point (after unit has stabilized)


The E02 and E01 errors should **NOT** be addressed by the end user. Switch the unit off and contact your VWR representative for repairs.

REPLACEMENT PARTS

See image on page 10.

DESCRIPTION	PART NUMBER
1. Lid assembly	880900-00
Lid	280632-00
Hinge	180108-00
Handle	180106-00
2. Lift	280516-00
3. Lift mount	180104-00
4. Locking acorn nut	180105-00
5. Fan guard	280620-00
6. Grill	280625-00
7. Lid gasket, front only	280605-00
8. Lid gasket, sides and back	280606-00
9. Safety switch	280509-00
10. Heater	380621-00
11. Thermostat	380723-00
12. RTD assembly	380622-00
13. Fan	280515-00
14. Insulation, left side	280600-00
15. Insulation, right side	280601-00
16. Rubber mat (33 x 27.9cm)	480004-00
17. Tray	280505-00
18. Knurled thumbscrew (#10-32 x 1.9cm)	180101-00
19. Motor	280633-00
20. V-belt	580000-00
21. Front panel overlay	380625-00
22. Power supply	380623-00
23. Display board	380624-00
24. Line filter	387022-00
25. IEC power entry module	386122-00
5 amp fuse (5 x 20mm)	380238-00
26. Feet	580002-00
Detachable 234cm power cord:	
120V	330100-00
EURO	330101-00
UK	330102-00
SWISS	330103-00

REPLACEMENT PARTS



ACCESSORIES

TEST TUBE RACKS - Requires Platform accessory.

DESCRIPTION	TUBE CAPACITY	PLATFORM CAPACITY	CAT. NO.
10 to 13mm Test Tube Rack	63	2	444-0278
14 to 16mm Test Tube Rack	48	2	444-0279
18 to 20mm Test Tube Rack	35	2	444-0280
22 to 25mm Test Tube Rack	24	2	444-0281
15mL Centrifuge Tube Rack	42	2	444-0283
50mL Centrifuge Tube Rack	12	2	444-0284
1.5 to 2mL Micro-Tube Rack	70	2	444-0276

Installation of Test Tube Racks:

1. Unit should be in standby mode.
2. Remove non-skid rubber mat. Slide Platform over-top included tray.
3. Using the hardware provided, attach brackets to the Platform. Position brackets as shown. The proper spacing between the brackets is 19.7cm.
4. Align the test tube rack side-to-side with the brackets on the Platform.
5. Press the test tube rack in place, you will hear the test tube rack clip into position. No mounting screws are required.

11

FLASK CLAMPS - Requires Platform accessory.

DESCRIPTION	MATERIAL	PLATFORM CAPACITY	CAT. NO.
10mL Erlenmeyer Flask Clamp	stainless steel	60	444-7040
25mL Erlenmeyer Flask Clamp	stainless steel	30	444-7041
50mL Erlenmeyer Flask Clamp	stainless steel	24	444-7042
125mL Erlenmeyer Flask Clamp	stainless steel	12	444-7043
250mL Erlenmeyer Flask Clamp	stainless steel	8	444-7044
500mL Erlenmeyer Flask Clamp	stainless steel	6	444-7045
1L Erlenmeyer Flask Clamp	stainless steel	4	444-7030
500mL Media Bottle Clamp	stainless steel	6	444-7035
1L Media Bottle Clamp	stainless steel	4	444-7036

Installation of Flask Clamps:

1. Unit should be in standby mode.
2. Remove non-skid rubber mat. Slide Platform over-top included tray.
3. Line up the hole(s) in the flask clamp to the hole(s) in the Platform. Using the screw(s) provided, hand tighten into place using a flathead screwdriver. **DO NOT** over tighten.
4. Insert flask/media bottle into clamp and wrap the spring around the neck of the clamp to hold secure.


ACCESSORIES

PLATFORM


Platform can be used for mounting Erlenmeyer flask/media bottle clamps and test tube racks. Slides over-top included tray. Tighten with adjustment screws. Stainless steel construction.

Description	Dimensions (L x W)	Part Number
Platform	33 x 27.9cm	444-7087


PLATFORM CONFIGURATIONS









Configurations of flask/media bottle clamps and test tube racks on the Platform accessory.

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