VWR Automatic vertical autoclaves

INSTRUCTION MANUAL

European Catalogue Number:

VAPOUR-Line^{eco} 50

481-0665

Version: 1.5
Issued: 2nd August 2016

Legal address of manufacturer

Europe

VWR International S.r.l.

Via San Giusto 85

20153 Milano (MI)

Italy

Tel.: 02-3320311/02-487791

<http://it.vwr.com>

Country of origin Italy

Table of contents

Use of the manual	4
Safety information	4
Precautions during installation	4
Unit description	5
Working principles	6
External component description	7
Control panel description	8
Operation	8
Sterilisation parameter settings	10
Configuration parameter settings	11
Sterilisation cycle	14
Validation	15
Maintenance and cleaning	16
Temperature/pressure chart	16
Troubleshooting	17
Electrical diagram	18
Printer	19
Accessories	20
Safety precautions	20
Warranty	21
Disposal	22

Use of the manual

Using this manual will help you to keep your instrument efficient. This manual should be kept in a place nearby the instrument for easy reference.

Safety information

The device must be correctly installed according to these instructions before starting any operation.

Never disconnect the device by pulling the electric cable from the socket.

Immediately replace the electric cable when it becomes damaged or worn.

When removing this item when connected to an electrical system pay particular attention to the connection cable.

Always disconnect from an electrical system when:

A. Repairing or maintenance; these operations must be carried out by qualified staff only

B. Cleaning the unit

Use original spare parts and accessories for any replacement.

Do not use this device in the presence of explosive gases.

Precautions during installation

When the electrical system socket is not compatible with device socket, replace it paying attention to the following coded cable colours:

Grey/Black/Brown: Live - Blue: Neutral - Yellow/green: Earth

Please ensure that the voltage and frequency of the electrical system is compatible with device.

The device must be connected to a ground wire. For any question please contact an electrician.

Before starting up the device ensure that the ventilation grids are free of any obstruction.

Unit description

General description

The VAPOUR-Line eco 50 autoclave operates by heating element and is classified as a “self-contained steam generator”.

Operating/project conditions

Working pressure 2,1 bar

Design pressure 2,2 bar

Hydraulic test pressure 3,3 bar

Hydraulic test temperature 20 °C

Overall dimensions

Height:	1005 mm	Height with open lid:	1300 mm
Width:	470 mm	Width with taps	550 mm
Depth:	550 mm	Depth with cable	570 mm
Weight:	60 kg		

Sterilisation vessel

Total nominal volume 50 l

Water required 4 l

Sterilisation vessel description

Vertical cylindrical chamber, manufactured in AISI 316L S/S electro polished.

AISI 316L lid, covered with insulating ABS to prevent burns, equipped with a handle for a friendly open/close operation.

Silicon cover gasket specially fitted for high temperatures.

Structure description

AISI 316L S/S square-shaped outer body.

Safety devices

The autoclave is equipped with the following safety devices:

1. Thermostatic over-heating cut-off
2. Microprocessor electronic safety; prevents total chamber overheating when temperature reaches 138 °C
3. Lid sensor
4. Electrically activated mechanical lock - prevents the lid opening if the temperature measured by the flexible or the vessel thermoprobe is higher than 80 °C. DO NOT FORCE THE DOOR OPEN. The door can be opened again only when both vessel and product temperature are below 80 °C
5. Shut-off safety valve when pressure gets over 2,2 bar

Hydraulic system

- Safety valve
- Wastewater pipe
- Pressure gauge
- Electrical heating element
- Draining tap for vessel water
- Hose connection for condensate discharge
- External port (rear panel) to allow the introduction of a certified thermoprobe
- Electro valve

Electrical system

Complies with CEI EN 61010-1, silicon fireproof cables.

- Heating element power: 3700 W
- Power: Triphase 380 V / 50 Hz
- Provided with tri-phase electric cable

Electric connections

Connections must be carefully followed:

Grey/Black/Brown - live

Blue - neutral

Yellow/green - ground wire

as per the enclosed wiring diagram

Installation requirements:

- Ground wired socket
- An industrial socket with a max. current up to 20 A with a residual current circuit breaker (RCCB) of 30 mA

Working principle

The sterilisation cycle consists of raising the chamber temperature to a desired value and keeping it for an established period of time.

It is possible to set the temperature value, bearing in mind that there is a correlation between temperature and pressure. Note that the autoclave is equipped with a pressure gauge which shows pressure measured in bar units.

The correlation between temperature and pressure in this autoclave is as follows:

121 °C - 1,03 bar

134 °C - 2,0 bar

Heating operation produce steam which helps air evacuation. The autoclave is equipped with an electro valve that automatically closes when the air contained in the vessel has been discharged.

When the total sterilisation cycle is completed an acoustic alarm will advise the operator; once the alarm rings the chamber starts its cool down operation.

External components description

Pos. 1 - Main switch

Pos. 2 - Vessel pressure gauge

Pos. 3 - Printer (standard)


Pos. 4 - Lid locking ring

Pos. 5 - Command panel

Pos. 6 - Lid opening handle


Pos. 7 - Steam drain hose nipple

Pos. 8 - Vessel water drain tap

Control panel description

The autoclave is equipped with specially designed microprocessor that controls the entire instrument and the sterilisation cycles.

Pos. 1 - Display

Multipurpose, it shows all the system messages and sterilisation cycle parameters such as temperature and time

Pos. 2 - Start button

Starts the sterilisation cycle shown on the display

Pos. 3 - Stop button

Stops the sterilisation cycle and resets all the alarms

WARNING! STOP button is an emergency button. Pushing it will abort the cycle and will open the valves. We don't recommend you push it unless a problem occurs during the cycle.

The same thing will happen if the instrument is turned off by the main switch.

Pos. 4 - Menu button

Enters the instrument setting menu

Pos. 5 - ▼ button

Scrolls up or increases the displayed value

Pos. 6 - ▲ button

Scrolls down or decreases the displayed value

Pos. 7 - Enter button

Confirms and stores the displayed value

Pos. 8 - LED

Reports the status of the autoclave

PRINTER ON: Lit when the printer is enabled


END CYCLE: Lit at the end of a sterilisation cycle

ALARM: lit when an alarm occurs

Operation

1. The autoclave should be placed in a room accessible by the operator only
2. Check that the heating element inside the vessel is fully covered by water
3. Check that the water tap (Pos. 8 page 7) is closed
4. Connect a suitable container to the hose connection of the condensate discharge port (Pos. 8 page 7)

5. Insert the basket with the items to be sterilised
6. Pull the safety ring placed on the left side of the cover to freely move the cover
7. Close the cover by using the handle
8. Turn on the autoclave by the main switch
9. The display will show the main menu:

10. Select the desired sterilisation cycle. If you want to modify the sterilisation parameter please follow the procedure described in the following chapters
11. Press START to start the cycle

Flexible thermo probe

The autoclave is equipped with a flexible thermo probe that allows temperature measurement in any point of

the chamber. The temperature measured by this probe is the PRODUCT TEMPERATURE. We suggest putting it in the coldest point in order to be sure that the sterilisation temperature has been reached. The probe can be put inside a small bottle containing liquid that has been sterilised. In such cases, please be aware that we recommend not to exceed 500 ml

The instrument has also a fixed thermoprobe that controls the temperature of the vessel. The temperature measured by this probe is the VESSEL TEMPERATURE.

You can stop the sterilisation cycle at any time by pressing the STOP button.

Stops the sterilisation cycle and resets all the alarms

WARNING! STOP button is an emergency button. Pushing it will abort the cycle and will open the valves. We don't recommend you push it unless a problem occurs during the cycle.

The same thing will happen if the instrument is turned off by the main switch.

Sterilisation parameter settings

1. Turn on the autoclave by the main switch (Pos. 1 page 7).
2. After the initial presentation the display will show Program 1 (P1)

P	1					1	2	1	C
							1	0	m

3. To modify this programme press the ENTER button, or press the ▲ or ▼ buttons until the programme to be modified is displayed and then press ENTER.

P	1				▶	1	2	1	C
							1	0	m

4. An arrow will be displayed next to the sterilisation temperature, press the ▲ or ▼ buttons until the desired sterilisation temperature is displayed and then press ENTER.

P	1					1	2	1	C
				▶			1	0	m

5. An arrow will be displayed next to the sterilisation time, press the ▲ or ▼ buttons until the desired sterilisation time is displayed and then press ENTER.
6. The modified data is automatically stored in the memory.
7. Press the START button to start the sterilisation cycle.

NOTES:

Program 5 has a double cycle making it possible to set two independent temperatures at the same time.

Program 6 can be used to run a cycle with the additional information of the F0 parameter. The calculation of this parameter is based on the setting selected in sub menu F0 in the option menu.

Program 7 can be used to run a cycle at the selected temperature until the selected value of the F0 parameter is reached. The calculation of this parameter is based on the setting selected in sub menu F0 in the option menu.


Configuration parameter settings

Press the MENU button to modify the instrument configuration parameter settings.


There are eight sub menus:

- **"Printer"** - modifies the printer settings
- **"F0"** - sets the value for F0 calculation on programmes 6 and 7
- **"Delay"** - sets a delayed start of the sterilisation cycle
- **"Print last cycle"** - prints the last sterilisation cycle completed
- **"Password"** - sets a password to prevent unauthorised modification of the sterilisation and instrument settings
- **"Date and time"** - modifies the autoclave date and time
- **"LCD contrast"** - modifies the adjustment of the display contrast
- **"Service"** - This menu is for maintenance authorised personnel only. The parameters stored in this menu can only be displayed but not modified...with the exception of the language.

Sub menu "Printer"

From this menu it is possible to enable or disable the printer plus set two different printing intervals for the sterilisation and heating/cooling phases. Press the ENTER button from the "PRINTER" menu; the workflow is as follows:

If the printer is enabled the “PRINTER ON” LED on the control panel will be lit.

Printing intervals can be set from a minimum of 1 minute up to a maximum of 30 minutes with a 1 minute step.


Sub menu “F0”

F0 parameter settings can be modified in this sub menu. The mathematical expression used for the F0 calculation (in minutes) is as follows:

$$F0 = \Delta t \cdot \sum 10^{\left(\frac{T - 121,11}{z}\right)}$$

where T is the vessel temperature, z is a microorganism-related parameter ranging from 4,0 to 15,0 (default 10), Δt is the time interval between measurements. The calculation of F0 starts after a threshold temperature T_s , ranging from 100 to 120 °C (default 100 °C), is reached.

From the F0 menu, press ENTER.


Values of z for some bacteria:

Microrganism	Mean value of Z parameter
Clostridium botulinum	110
Bacillus stearothermophilus	6
Bacillus subtilis	10
Bacillus megaterium	7
Bacillus cereus	10
Clostridium sporogenes	13
Clostridium histolyticum	10

Sub menu "Delay"

▲	D	E	L	A	Y
▼	0	0	:	0	0

From this menu it is possible to set a delayed start for the sterilisation cycle. From the "DELAY" screen press the ENTER button:

Once the delay has been set, after pressing the START button, the display will show the delay and the remaining time before the cycle begins. The delay is stored in the memory, set 00:00 as delay to reset the delay. The maximum settable delay is 23:59

Sub menu "Print last cycle"


▲	P	R	I	N	T	L	A	S	T
▼	C	Y	C	L	E				

From this menu it is possible to print the last completed cycle by pressing the ENTER button.

Sub menu "Password"

▲	P	A	S	S	W	O	R	D
▼								

From this menu it is possible to set a password to prevent unauthorised modification of the sterilisation and instrument settings. The password can be set by pressing the ENTER button from "PASSWORD".

From the main menu the password will be requested every time the ENTER or MENU buttons are pressed. A password will not be requested if the password value has been set to 0000 (default value).

Sub menu "Date and time"

▲	D	A	T	E	A	N	D	T	I	M	E
▼											

From this menu it is possible to adjust the date and time. Press the ENTER button from the "DATE AND TIME" menu.

02/03/2010	▲ o ▼ to increase or decrease
00:00	ENTER to skip to the next field
	MENU to exit

Sub menu "LCD contrast"

▲	L	C	D	C	O	N	T	R	A	S	T
▼				5	5	[	%	]			


From this menu it is possible to adjust the display contrast. Press the ENTER button from the menu "LCD CONTRAST".

The contrast default value is set to 55%.

▲ o ▼ to increase or decrease
MENU to exit


Sub menu "Service"

This menu is for maintenance authorised personnel only. The parameters stored in this menu can only be displayed but not modified... with the exception of the language.

Sterilisation cycle

1. To start a cycle, select the desired sterilisation programme
2. Press the START button
3. The display will show on the upper row the selected cycle (temperature and time), while on the lower row will show the actual product temperature and the remaining sterilisation time. In order to see the vessel temperature, press MENU and use the arrows in order to see both Vessel and Product temperature.

4. Once the sterilisation temperature has been reached, the sterilisation time count will begin. At the end of the sterilisation time a buzzer is activated
5. In the cooling process, when the vessel temperature goes below 100 °C the LED END CYCLE will light up
6. Once the cooling process is completed, open the lid with the handle
7. Switch off the autoclave via the main switch (Pos. 1 page 7)

NOTE:

Press the STOP button to end a running cycle.

WARNING! STOP button is an emergency button. Pushing it will abort the cycle and will open the valves. We don't recommend you push it unless a problem occurs during the cycle.

The same thing will happen if the instrument is turned off by the main switch.

The set values are stored in the memory will not be lost after switching the instrument off.

Validation

The autoclave should be checked each time with chemical indicator and validated monthly with biological indicators.

Maintenance and cleaning

The autoclave should be inspected annually according to "Good Laboratory Practice".

Daily cleaning of the autoclave will prevent most future faults.

Clean the autoclave with a mild detergent using a wet soft cloth. Replace the vessel water daily with new, fresh water. Long-term use of dirty water can cause premature failure of the heating element.

WARNING

Do not use any cleaning product containing halogens or other aggressive substances to clean the autoclave as this can cause rust.

Temperature/pressure chart

Temperature (°C)	100	102	105	109	112	115	119	120	128	134
Pressure (bar)	0	0,1	0,2	0,4	0,5	0,7	0,9	1,0	1,5	2,0

Troubleshooting


A malfunction or an alarm condition can happen during normal working procedures.

In the following table the most common fault conditions are described as well as their possible corrective actions.

BEFORE OPENING THE DOOR, PLEASE CHECK THE PRESSURE GAUGE THAT THERE IS NO PRESSURE INSIDE THE VESSEL. IN SUCH A CASE LET THE CHAMBER COOL DOWN IN ORDER TO RELEASE THE EXCESS PRESSURE. DO NOT FORCE THE LID OPEN!!!

EFFECT	CAUSE	ACTION
Main switch and display are off	No power to the autoclave	- Check mains voltage - Check power connections - Check fuses
Display is off and the main switch LED is on	No power to the main board	Contact technical services, possible faulty main board or low voltage power supply
Lid does not Open	Product temperature is higher than 80 °C	Wait until product temperature is below 80 °C. DO NOT FORCE THE LID OPEN
Lid does not Open and the product temperature on the display is below 80 °C	Vessel temperature is higher than 80 °C. Press MENU and arrow to check it.	Wait until vessel temperature is below 80 °C. DO NOT FORCE THE LID OPEN
Heating element is not working	- Faulty heating element - Vessel overheating, the safety thermostat is cutting off the heater	Switch off the autoclave, wait for ten minutes and then turn on it back on; if the problem persists contact technical services
The autoclave is not reaching the programmed temperature	Faulty electro valve	If the problem persists contact technical services
The autoclave only reaches 100/110 °C, but not over this temperature	Faulty electro valve	Contact technical services
The vessel temperature is more than 135 °C	Faulty solid state relay	Contact technical services
ALARM LED is lit and the display shows: ALARM DOOR OPEN	Door sensor is not activated	Fully move the closing handle to its closed position; if the problem persists contact technical services
ALARM LED lit and the display shows: ALARM OVERTEMPERATURE	Vessel overheating	Contact technical services
Display shows "alarm CRC"	No water inside the chamber	Add water
Display shows "alarm CRC", but there's water inside the chamber	In order to work properly the water sensor, as it's based on conductivity, must be in contact with water with a minimum amount of salt	Add a small spoon of salt or a glass of tap water to the autoclave

Electrical diagram

- | | | |
|--------------------------|---------------------------|---------------------------|
| Pos. 1 Line fuses | Pos. 6 Electronic board | Pos. 11 Solid state relay |
| Pos. 2 Main switch | Pos. 7 Lid micro switch | Pos. 12 Heating element |
| Pos. 3 Noise filter | Pos. 8 Electro valve | Pos. 13 Safety relay |
| Pos. 4 24 V power supply | Pos. 9 Printer | |
| Pos. 5 Thermo probe | Pos. 10 Safety thermostat | |


Printer

The autoclave is equipped with a connection to a printer.

The stored information is as follows:



- a) Autoclave model and serial number
- b) Date and time of the sterilisation cycle beginning
- c) Sterilisation parameters set (temperature and time)
- c) Vessel temperatures recorded with the set interval time

It is also possible to write down the batch number and operator sign.

- Pos. 1 Opening button
- Pos. 2 Paper roll container
- Pos. 3 Printing head
- Pos. 4 Status LED
- Pos. 5 Paper feed button



How to insert a new paper roll

To open the printer, press the central green button.

Place a new paper roll according to the procedure described in the following pictures.

If the paper has been correctly placed, the LED in the central green button will remain steady. A flashing LED means an error. See paragraph 8.1 of this manual for the directions to modify the printing interval.

Accessories

Description	Pk	Cat. No.
ECObox-50- basket 320x285	2	481-0672
ECObox-50 mini basket 310x190	3	481-0668

Safety precautions

Opening and closing

The autoclave has a safety block when the lid is open. This block will prevent the lid falling on the hands of the operator. ALWAYS make sure that the block is working properly and that there's sufficient space between the instrument and the wall behind, in order to allow the lid to open completely and block to be activated.

After a cycle, the instrument should only be opened when temperature is below 80 °C and the inner pressure is at 0. Autoclaves have safety devices in order to avoid accidental opening. In this case DO NOT FORCE the instrument open.

Opening the autoclave when power is off is NOT recommended.

Personal Safety Systems (PSS)

All the safety devices on this autoclave will prevent accidental opening. However, sterilised loads can undergo thermal shock during a cycle that can cause breakages during opening, although in very rare circumstances.

The following PSS workwear should be near the instrument and always worn by personnel:

- Heat resistant gloves
- Heat resistant face mask (e.g. made of polycarbonate)
- Protective glasses
- Labcoat

Sterilisation of liquids

Please note that liquids in a load have a thermal capacity bigger than the steam in the chamber. For this reason they're always a "step behind" the temperature measured in the chamber. The more viscous or bulky the liquid is, the more this phenomenon may occur.

So, when the chamber is less than 100 °C, liquid inside the load can be at a temperature higher than this.

Consequently if the autoclave is opened, liquid inside (that is at a much higher temperature) may come out in a turbulent manner and dirty the internal chamber or, (worst case) splash the operator. If the container is closed it can explode.

If you're sterilising liquids, always consider **at least** a 20 °C difference between chamber and temperature and load temperature. The instrument should therefore be opened at 50 - 60 °C.

We strongly recommend **NOT** to sterilise liquids in containers bigger than 2 litres. Containers should be no more than $\frac{3}{4}$ full.

Warranty

VWR International warrants that this product will be free from defects in material and workmanship for a period of two (2) years from date of delivery. If a defect is present, VWR will, at its option and cost, repair, replace, or refund the purchase price of this product to the customer, provided it is returned during the warranty period. This warranty does not apply if the product has been damaged by accident, abuse, misuse, or misapplication, or from ordinary wear and tear. If the required maintenance and inspection services are not performed according to the manuals and any local regulations, such warranty turns invalid, except to the extent, the defect of the product is not due to such non performance.

Items being returned must be insured by the customer against possible damage or loss. This warranty shall be limited to the aforementioned remedies. IT IS EXPRESSLY AGREED THAT THIS WARRANTY WILL BE IN LIEU OF ALL WARRANTIES OF FITNESS AND IN LIEU OF THE WARRANTY OF MERCHANTABILITY.

Compliance with local laws and regulations

The customer is responsible for applying for and obtaining the necessary regulatory approvals or other authorisations necessary to run or use the product in its local environment. VWR will not be held liable for any related omission or for not obtaining the required approval or authorisation, unless any refusal is due to a defect of the product.

Equipment disposal

This equipment is marked with the crossed out wheeled bin symbol to indicate that this equipment must not be disposed of with unsorted waste.

Instead it's your responsibility to correctly dispose of your equipment at lifecycle end by handing it over to an authorised facility for separate collection and recycling. It's also your responsibility to decontaminate the equipment in case of biological, chemical and/or radiological contamination, so as to protect from health hazards the persons involved in the disposal and recycling of the equipment.

For more information about where you can drop off your waste equipment, please contact your local dealer from whom you originally purchased this equipment.

By doing so, you will help to conserve natural and environmental resources and you will ensure that your equipment is recycled in a manner that protects human health.

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VWR International AB
VWR International AG
Teknolojileri Ltd.Şti.
VWR International Ltd