VWR[®] Microbiological Air Sampler SAS Super IAQ

INSTRUCTION MANUAL

European catalogue numbers:

710-0868

710-0919

Version: 1.2
Issued: July 2016

Legal address of manufacturer

Europe

VWR International Srl

Via S. Giusto 85

20153 Milano

Italy

Tel.: + 39 02-3320311

<http://it.vwr.com>

Country of origin: Italy

Table of contents

Legal address of manufacturer.....	2
Warning.....	5
Safety information.....	5
Kit and unit description	6
Technical features	7
Intended use	7
Principle	7
The basic idea.....	7
Brief operating instructions	8
Brief operating instruction for SAS Super IAQ.....	8
Installation	8
Practical use of contact plates	8
List of menus and utility sub menus	9
Preliminary inspection.....	9
Holder adjustment.....	9
Petri head adapter (optional)	9
90 mm Petri dish filling.....	10
Tripod installation (optional).....	10
Functions	10
Software	10
Starting with the same air volume as the previous sample	10
“USER MODE”	10
“PROGRAM VOLUME” (modification of the value of a stored volume)	11
“DELAY MODE”	12
“UTILITY MODE”	12
“SET TIME”	13
“SET AUTOSWITCH”	13
“DISPLAY RECORD”	14
“CLEAR RECORDS”	14
“LANGUAGE”	14

“IDENTIFY”	15
“SAMPLING SITE”	15
List of system messages	15
Accessories	16
Troubleshooting	17
Technical service	18
Warranty.....	18
Compliance with local laws and regulations	18
Equipment disposal.....	19

Warning

This document is the property of VWR International PBI S.r.l. - Milan - Italy

It may not be duplicated or distributed without the owner's authorisation.

Patent pending

Manufactured by VWR International PBI S.r.l. - Milan, Italy

Safety information

Use the device only for the indicated purposes.

The device must be correctly used according to this instruction manual before starting any operation.

Immediately replace any electrical cable when damaged.

Always disconnect the charger before:

- Repairing or maintenance; these operations must be carried out by qualified staff
- Cleaning the unit

Use original spare parts and accessories for any replacement.

Do not use this device in the presence of explosive gas.

Please follow the guidelines below and read this manual in its entirety to ensure safe operation of the unit.

Be aware that the voltage and frequency of the electrical system are compatible with the electrical requirements of the battery charger.

Never use a non OEM charger to charge the air sampler.
Use of improper charger may damage the unit.

Kit and unit description

1. Soft carrying case
2. SAS IAQ control unit
3. Aluminium head
4. Battery charger
5. Remote control
6. "Clear" button
7. "Enter" button
8. Arrow "Up"

9. "Start" button
10. Arrow "Down"
11. Standing feet
12. Printer connection (not used)
13. Main switch
14. Serial number
15. CE mark
16. Charger connection

Technical features



Model	SAS IAQ Contact	SAS IAQ Petri
Cat. No.	710-0868	710-0919
Airflow	100 l/min	100 l/min
Use with	55 mm contact plates	90 mm Petri dishes
Portable	Yes	Yes
Battery	Rechargeable	Rechargeable
Battery life	70 000 l	70 000 l

Intended use

Principle

The Surface Air System (SAS) encompasses several models which use the same principle. Air is aspirated at a fixed speed for variable time through a head which has been machined with a series of specially designed small holes. The resulting laminar airflow is directed onto the agar surface of a “RODAC plate” (or Petri dish) containing medium for microbiological analysis.

When the selected sampling cycle is completed, the plate is removed and incubated. The organisms are then visible to the naked eye and can be counted for the assessment of contamination levels.

The basic idea

The key features of the Surface Air System (SAS) are:

1. To use a simple and inexpensive “contact plate” (“SURFAIR PLATE”, “RODAC”) for surface, hand or air control. These plates are very well known and easily available. They can be purchased ready poured with different media.
2. To sample a known volume of air for a variable time to provide a range of sampling volumes.
3. To aspirate air in a laminar flow pattern with sufficient velocity to impact organisms onto an agar surface.
4. To accumulate data on the level of hygiene in each environment so that fluctuations can be monitored.
5. To take advantage of advanced electronics for more reliable results in different operating conditions.
6. To have the flexibility to choose between 55 mm contact plates or 90 mm standard Petri dishes.


Brief operating instructions

Brief operating instruction for SAS Super IAQ

1. Press main switch.
2. Press 'START' button to sample the same air volume as the last sampling cycle.
3. To change the air volume, use the arrow buttons and select 'USER MODE' from the menu.
4. Refer to the instruction manual to change settings.
5. Press 'ENTER' to confirm selection.
6. Press arrows for menu selection to modify other pre-set parameters.
7. Press 'ENTER' to reach the sub menu of the chosen parameter.
8. Refer to the instruction manual to change settings.
9. Press 'ON/OFF' switch at the end of sampling operations.
10. Press 'CLEAR' each time you need to end an action. The unit will then return to the initial configuration.

Installation

Practical use of contact plates

1. Remove the aspirating head.
2. Insert an identified, closed and prepared plate and remove the plate lid.
3. Replace the aspirating head.
4. Select required air volume and start the unit. The airflow is directed onto the agar surface of the plate.
5. At the end of the cycle, remove the aspirating head.
6. Close and remove the plate.
7. Incubate.
8. Count the colonies, record the results on the microbiological air sampling report and read the results.

List of menus and utility sub menus

START FOR XXX	This message indicates the air sampler is ready to sample the volume shown (XXX) after pressing the 'Start' button
USER MODE	Choose from 8 programmable sampling volumes
PROGRAM MODE	Modify the 8 programmable sampling volumes
DELAY MODE	Set a delayed start
UTILITY MODE	This leads to the following sub menu:
	SET TIME To adjust the date and time
	SAMPLING SITE To identify the sampling point
	IDENTIFY To identify the operator's name
	LANGUAGE To display text in different languages
	CLEAR RECORD To clear the recorded sampling data
	DISPLAY RECORD To show the recorded sampling data
	AUTO SWITCH OFF To disconnect the automatic switch-off when the "Infrared Remote" is used

Preliminary inspection

The instrument is subject to specific quality tests before shipment and is carefully packed to avoid possible damage during transit. However, a visible check should be carried out as soon as possible to determine any damage. This must be reported immediately. The following procedures should be followed to check that the unit is working properly.

The instrument battery pack must be charged for at least 4 hours before running the test.

After the 'ON/OFF' button is pressed an automatic visual presentation will be displayed.

Holder adjustment

Contact plate and Petri dish holders can be adjusted (using a screw driver) in case the available plate is slightly different in diameter from the standard 55 mm contact plate or 90 mm Petri dish.

Petri head adapter (optional)

An adapter is available in aluminium or stainless steel which allows the use of standard 90 mm Petri dishes with SAS for contact plates.

90 mm Petri dish filling

We recommend you fill the standard 90 mm disposable Petri dish with no more than 18 - 20 ml of agar to avoid the medium touching the inside surface of the aspirating head.

Tripod installation (optional)

The SAS SUPER IAQ can be fixed to a table, floor tripod or wall support. The screw thread connection is located under the unit, between the two front feet.

Functions

Software

Press the 'ON/OFF' switch (black button) to switch the unit on.

The 'ON/OFF' switch (black button) must be pushed twice to switch on the unit if the automatic switch-off timer was used previously.

Press 'CLEAR' at any time to return to the initial 'START FOR' main screen.

The air sampler is provided with eight programmable air volumes ('User Mode') selectable by the operator.

Volume selection is made by pressing the "▲" or "▼" arrows when the programme is in the relevant sub menu.

Suggested air volumes:

- Contaminated areas (communities, processing rooms, etc) 10 - 200 litres of air
- Normal areas (laboratory benches, houses, etc) 200 - 500 litres of air
- Sterile or high risk areas (cleanrooms, operating theatres, etc) 500 - 1000 litres of air

Starting with the same air volume as the previous sample

Switch on the unit by pushing the main switch. You will then see the following display with the last aspirated volume of air:

START FOR XXX
MM/DD/YY - HH:MM

Push the 'START' button to start a sampling cycle with the displayed volume of air.

[XX] → XXX
MM/DD/YY - HH:MM

'USER MODE'

Up to eight additional volumes can be programmed into the unit. User-selectable volumes can be selected from the 'USER MODE', while the selectable volumes can be modified using the PROGRAM MODE' function described in the next paragraph.

Switch on the instrument and wait until the main screen is displayed:

START FOR XXX
MM/DD/YY - HH:MM

Press the “Up” or “Down” arrows to reach the:

USER MODE
MM/DD/YY - HH:MM

Press the ‘ENTER’ button, now you can scroll the eight user-defined sampling volumes.

S.PROG XXXX
MM/DD/YY - HH:MM

Once the desired sampling volume has been reached press the ‘ENTER’ button to confirm.

Press the ‘START’ button to begin the sampling cycle:

START FOR XXX
MM/DD/YY - HH:MM

There are a total of eight available “USER MODE’ programmes.

The maximum volume of air for each sampling cycle is 1999 litres.

‘PROGRAM VOLUME’ (modification of the value of a stored volume)

Using this procedure, it is possible to memorise up to eight different volumes (from 1 to 1999 litres of air). Switch on the instrument and wait until the main screen is displayed.

Press the “Up” or “Down” arrows to reach the:

START FOR XXX
MM/DD/YY - HH:MM

Press the ‘ENTER’ button, now you can scroll the eight user-selectable sampling volumes.

Select a volume by pressing the “UP” or “DOWN” arrows and then pressing ‘ENTER’:

PROGRAM MODE
MM/DD/YY - HH:MM

Once the sampling volume has been set, press the 'ENTER' button to confirm.

S.PROG XXXX
MM/DD/YY - HH:MM

The display will show the old and new volume.

Press the "Up or "Down" arrows to select the first digit and press 'ENTER' to confirm.

XXX → 0000
MM/DD/YY - HH:MM

Press the Up or "Down" arrows to select the second digit and press 'ENTER' to confirm.

XXX → X000
MM/DD/YY - HH:MM

'DELAY MODE'

With this function it is possible to delay the beginning of the sampling cycle. If a delay has been selected, the sampling cycle will begin after the selected delay time.

Switch on the instrument and wait until the main screen is displayed:

START FOR XXX
MM/DD/YY - HH:MM

Press the "Up" or "Down" arrows to reach the:

DELAY MODE
MM/DD/YY - HH:MM

Press the 'ENTER' button and the delay time is displayed:

DELAY 01 MIN
MM/DD/YY - HH:MM

Select the desired delay time using the "Up or "Down" arrows. The selectable delay values are 1, 2, 3, 5, 10 or 20 minutes. Press the 'ENTER' button to confirm.

When a delay has been selected, after the "START" button has been pressed, a delay warning message will be displayed.

** D E L A Y **

'UTILITY MODE'

With this function it is possible to modify all the settings of the air sampler.

The 'UTILITY MODE' is divided into ten sub menus:

Set time; Set Autoswitch; Display Record; Clear Record; Language; Identify; Sampling site.

Always start from 'UTILITY MODE' to reach one of these sub menus.

Switch on the instrument and wait until the main screen is displayed:

START FOR XXX
MM/DD/YY - HH:MM

Press the "Up" or "Down" arrows to reach the:

UTILITY MODE
MM/DD/YY - HH:MM

Press the "ENTER" button, the first sub menu is then listed.

'SET TIME'

This option is used to programme day, month, year and time of day.

SET TIME
MM/DD/YY - HH:MM

Press 'ENTER' and the "Up" or "Down" arrows to change the month

Press 'ENTER' and the "Up" or "Down" arrows to change the day

Press 'ENTER' and the "Up" or "Down" arrows to change the year

Press 'ENTER' and the "Up" or "Down" arrows to change the hours

Press 'ENTER' and the "Up" or "Down" arrows to change the minutes

Press 'ENTER' to confirm and exit

'SET AUTOSWITCH'

To save battery consumption the SAS automatically switches off after 4 minutes. The 'SET AUTOSWITCH' option is used to enable or disable the automatic switch-off.

From the 'UTILITY MODE' select 'SET AUTOSWITCH':

SET AUTOSWITCH
MM/DD/YY - HH:MM

Press 'ENTER' and the "Up" or "Down" arrows to change between the two options:

AUTOSWITCH ON
MM/DD/YY - HH:MM

AUTOSWITCH OFF
MM/DD/YY - HH:MM

Select the desired options and press 'ENTER' to confirm and exit.

This setting will be lost after the instrument is switched off.

‘DISPLAY RECORD’

The last 99 samples are memorised in the file ‘DISPLAY RECORD’. Each sample is identified in chronological order and shows the date, time, operator, site and volume of air sampled.

From the UTILITY MODE select ‘DISPLAY RECORD’.

DISPLAY RECORD
MM/DD/YY - HH:MM

Press ‘ENTER’ to see the memorised values:

001 MM/DD HH:MM
ID. /SITE VOL. OP

The following parameters are recorded:

Progressive number; Month and Day; Time; Identification/Site; Sampling Volume.

If there is no data stored in memory the display will show:

‘CLEAR RECORDS’

This option is used to delete all the data memorised in the DISPLAY RECORD. Before starting this procedure, please be certain that existing data is not required or that it has been downloaded.

From the UTILITY MODE select ‘CLEAR RECORD’.

CLEAR RECORD
MM/DD/YY - HH:MM

Press ‘ENTER’ to delete all the stored data, the display will show the resetting status:

RESETTING...
▶▶▶▶

‘LANGUAGE’

Menu text may be selected by choosing from different languages.

From the UTILITY MODE select ‘LANGUAGE’:

LANGUAGE
MM/DD/YY - HH:MM

Press 'ENTER' to list all the available languages:

ENGLISH

Use the "Up" or "Down" arrows to select the desired language and press 'ENTER' to confirm.

There are six available languages:

English, French, Spanish, German, Portuguese and Italian.

'IDENTIFY'

This option is used to identify the operator. This should be changed if different operators use the sampler and especially if the data is to be printed.

From the UTILITY MODE select 'IDENTIFY'.

SAMPLING SITE
MM/DD/YY - HH:MM

Press 'ENTER' to modify identification code:

IDENTIFY
XXXX

Select the desired character using the "Up" or "Down" arrows and confirm with 'ENTER'.

Select the next desired characters and confirm them by pressing 'ENTER'.

The selected values are stored in the memory and will be kept even after the instrument has been switched off.

'SAMPLING SITE'

Site identification should be changed for samples taken at different sites, especially if the results are to be printed.

From the UTILITY MODE select 'SAMPLING SITE'.

SAMPLING SITE
MM/DD/YY - HH:MM

Press 'ENTER' to modify:

SAMPLING SITE
XXXX

Select the desired character using the "Up" or "Down" arrows and confirm with 'ENTER'. Select the next desired characters and confirm them by pressing 'ENTER'.

The selected values are stored in the memory and will be kept even after the instrument has been switched off.

List of system messages

Besides the messages described in the previous paragraph, the following messages can be displayed:

LOW BATTERY	Battery must be recharged
MOTOR ERROR	Motor not working or not connected
- OFF	The instrument is switching off

Accessories

Description	Cat. No.
Spare parts	
Aluminium aspirating head for contact plates, Ø 55 mm	710-0892
Aluminium aspirating head for Petri dishes, Ø 90 mm	710-0886
SAS Super universal battery charger	710-0993
Remote control	710-0962
Soft carrying case	710-0896
Other aspirating heads for contact plates, Ø 55 mm	
Stainless steel aspirating head for contact plates, Ø 55 mm	710-0880
Sterile daily head for contact plates, Ø 55 mm	710-0890
Other aspirating heads for Petri dishes, Ø 90 mm	
Stainless steel aspirating head for Petri dishes, Ø 90 mm	710-0878
Sterile daily head for Petri dishes, Ø 90 mm	710-0891
Other accessories	
Aluminium carrying case	710-0875
Bio-Transport autoclavable carrying case	113-8185
Handle for Bio-Transport carrying case	113-8186
Floor tripod	710-0889
SAS-Holder table and wall stainless steel	710-0963
Adaptor* to convert contact plate model to accept 90 mm Petri dishes	710-0882
SAS stainless steel Petri head + adaptor	710-0877
SAS aluminium Petri head + adaptor	710-0879

Troubleshooting

Review the information in the table below to troubleshoot operating problems.

Problem	Cause	Solution
Unit does not start	Battery is low	Recharge battery
	Battery charger is not working	Check battery charger and replace if required (Cat. No. 710-0993)
	Battery is too old	Check battery and replace if required
Battery discharges after few minutes working	Battery is low	Recharge battery
	Battery charger is not working	Check battery charger and replace if required (Cat. No. 710-0973)
	Battery is too old	Check battery and eventually replace it
“LOW BATTERY” message	Battery is low	Recharge battery
Infrared remote switch doesn't switch on the unit	Remote control battery is low	Replace infrared remote battery
Microbiological media is dehydrated after sampling	Media is corrupted	Check expiration date of media and if the agar was not dehydrated before sampling
	Sampling time is too long	Shorten the sampling time

Technical service

Web resources

Visit the VWR website at vwr.com for:

- Complete technical service contact information
- Access to the VWR online catalogue and information about accessories and related products
- Additional product information and special offers

Contact us for information or technical assistance contact your local VWR representative or visit vwr.com.

Warranty

VWR International warrants that this product will be free from defects in material and workmanship for a period of two (2) years from date of delivery. If a defect is present, VWR will, at its option and cost, repair, replace, or refund the purchase price of this product to the customer, provided it is returned during the warranty period. This warranty does not apply if the product has been damaged by accident, abuse, misuse, or misapplication, or from ordinary wear and tear. If the required maintenance and inspection services are not performed according to the manuals and any local regulations, such warranty turns invalid, except to the extent, the defect of the product is not due to non performance.

Items being returned must be insured by the customer against possible damage or loss. This warranty shall be limited to the aforementioned remedies. IT IS EXPRESSLY AGREED THAT THIS WARRANTY WILL BE IN LIEU OF ALL WARRANTIES OF FITNESS AND IN LIEU OF THE WARRANTY OF MERCHANTABILITY.

Compliance with local laws and regulations

The customer is responsible for applying for and obtaining the necessary regulatory approvals or other authorisations necessary to run or use the product in its local environment. VWR will not be held liable for any related omission or for not obtaining the required approval or authorisation, unless any refusal is due to a defect of the product.

Equipment disposal

This equipment is marked with the crossed out wheeled bin symbol to indicate that this equipment must not be disposed of with unsorted waste.

Instead it's your responsibility to correctly dispose of your equipment at lifecycle end by handing it over to an authorised facility for separate collection and recycling. It's also your responsibility to decontaminate the equipment in case of biological, chemical and/or radiological contamination, so as to protect persons involved in the disposal and recycling of the equipment from health hazards.

For more information about where you can drop off your equipment waste, please contact your local dealer from whom you originally purchased this equipment.

By doing so, you will help to conserve natural and environmental resources and you will ensure that your equipment is recycled in a manner that protects human health.

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VWR International AB
VWR International AG
Teknolojileri Ltd.Şti.
VWR International Ltd