Quick Start Guide

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CONTENTS

1	Safety instructions	2
2	License notice	2
3	System overview	3
4	Installation	3
4.1	Content of delivery	3
4.2	Installing the thermal cycler	3
4.3	Operating the thermal cycler	3
5	Cycler operation	4
5.1	Start screen = 'Run'	4
5.2	Running stored programs ('Run'):	4
5.3	Incubation mode ('Run')	5
5.4	Create or edit PCR programs ('Programs')	5
5.5	Further dialogues: Diagnostics, GLPs, System	9
5.6	Further questions...?	10
6	Technical service and ordering information	11
7	Technical data	12
8	Maintenance and repair	14
8.1	Cleaning, maintenance and repair	14
8.2	Over temperature test	14
8.3	Temperature check-up	14
8.4	Calibration	14

A more detailed instruction manual is available
for download on our website vwr.com.

1 SAFETY INSTRUCTIONS

Before the first use of the thermal cycler please read the entire instruction manual.

Special notice should be taken of the following:

Caution - dangerous voltage

Please ensure the voltage indicated on the device exactly matches your local electrical supply.

Caution - dangerous explosive material

Explosive or reactive material mustn't be heated in the thermal cycler.

Caution - liquids

Ensure that no liquids can enter the device.

Caution - hot surfaces

The thermal plate and the heated lid quickly reach temperatures above 50 °C. There is danger of burning! Keep the lid closed until the temperatures reach 30 °C or less.

Only use materials (tubes, tube caps, plates and sealing films) which are heat resistant to 120 °C.

Caution – environment

The devices ventilation system must not be covered.

Danger – line voltage

Incorrectly grounded products can be extremely hazardous in the event of a fault.

Use only a 3-conductor power cable (3x 0.75 mm²) with protective ground. The power connector may only be plugged into a socket with a protective ground. The protection must not be nullified by an extension cable without protective ground.

2 LICENSE NOTICE

This instrument is licensed for research and development and for uses other than human *in vitro* diagnostics under one or more of the following patents of Applied Biosystems Corporation: U.S. Patent Nos. 5,656,493, 5,038,852, 5,333,675, 5,475,610 (claims 1-159 and 164-166), 6,703,236 (claims 7-10) and 7,238,517 or corresponding claims in their non U.S. counterparts. No right is conveyed expressly, by implication or by estoppel under any other patent of Applied Biosystems, including but not limited to U.S. Patent No. 6,814,934 and its non U.S. counterparts, which describe and claim thermal cyclers capable of real time detection.

3 SYSTEM OVERVIEW

VWR thermal cyclers are the ideal PCR devices for all molecular biological and biochemical laboratories working in the field of basic research or routine diagnostic. The 250 °C HTR (High Temperature Range) Peltier technology allows quick changes in temperature up to 5 °C/sec (max). The individual monitoring and controlling of each of the eight or 16 Peltier elements per block equipped with 'Long Life Technology' ensure the tightest correlation of target and actual temperature across the entire thermal plate leading to highest reproducibility of the results. To avoid condensation each block is equipped with a heated lid.

The VWR thermal cyclers are equipped with a touch-sensitive, graphic coloured TFT display of the newest generation, allowing operation of the device by simply touching symbols and graphic elements on the screen. This leads to easy and intuitive handling and programming of the device. Additionally a mouse can be connected via USB (reboot required). Furthermore an Ethernet port is available to connect the devices with Microsoft Windows® or Linux networks. It is therefore possible to use an external server for archiving programs and GLP reports which will be available globally in the network (e.g. for other thermal cyclers). Transfer of data can be done by using commercially available USB memory sticks enabling a nearly unlimited saving space for programs and GLP reports, which can be printed using a printer connected to the thermal cycler via the network.

Note: The USB ports support only standard cable mice and standard USB sticks.
The maximum acceptable connecting cable length of the port is 3 m!

4 INSTALLATION

4.1 Content of delivery

- 1x Thermal cycler
- 1x Power cable
- 1x Short manual

(A detailed instruction manual in English or German is available for download on our website vwr.com)

- 1 CD PC software: VWR PCR Cycler Master software

Please check delivery for completeness and transport damage upon arrival. If any transport damage is noticed, contact your distributor or manufacturer immediately.

NOTICE: IF THERE IS ANY DAMAGE TO THE DEVICE DO NOT USE!

4.2 Installing the thermal cycler

Check the packaging for any transport damage. Remove the packaging and place the thermal cycler onto a solid surface. The device should not be exposed to direct sunlight. There should be enough space to make sure that the ventilation grills on the front and back are not covered and sufficient ventilation is provided. The accessibility of the power switch must not be hindered. There should be a distance of at least 25 cm (10 inch) to the wall or the next device. Two thermal cyclers should not be placed back-to-back or back to front. Room temperature between 4 °C / 39 °F and 25 °C / 77°F is advised.

All original packaging should be stored in case a need arises to return the instrument to the supplier.

4.3 Operating the thermal cycler

Before starting up the thermal cycler for the first time make sure that the voltage range of the device is suitable to the electricity supply. Plug in the mains plug of the thermal cycler to the electricity supply. The outlet and the power switch are located at the back of the device. The serial number is written onto a typed label at the bottom of the device.

The device should only be used by authorised staff. To ensure this, a user administration option with appointing rights for single users is available in the device.


5 CYCLER OPERATION

The operation of the thermal cycler can be run by the internal touch screen, either using the touch-sensitive surface of the display or by using a mouse connected to a free USB port.

If the user has to enter numbers or text the corresponding field has to be touched on the touch screen or clicked with the mouse. An appropriate keyboard will appear on the screen according to the expected entry.

5.1 Start screen = 'Run'



The last PCR program run is displayed in the 'Run' dialogue and can directly be restarted with the 'Start PCR' button.

Note: With the VWR Doppio double block thermal cyclers the desired block can be chosen by selecting 'L' (left) and 'R' (right). This option is not available for the VWR UNO⁹⁶, UNO³⁸⁴ and Ristretto.


5.2 Running stored programs ('Run'):

1. Open the selection list in the 'Run' option using the  button.
2. Select the desired program in the selection list and load the 'Run' option with the 'Open' button.
3. Start the loaded program by using the  button.



5.3 Incubation mode ('Run')

Use the  button in the 'Run' option to set the device to incubation mode. The temperature for the thermal plate and lid can be set manually. The clock shows the time elapsed in the incubation mode.

5.4 Create or edit PCR programs ('Programs')

1. To create or edit PCR programs select the 'Programs' option.
2. By touching the 'Programs' option for the first time during a session, the selection dialogue shows where a new program can be selected ('New') or an existing program can be opened ('Open') for editing.
3. On the left side of the program editor all available commands are displayed. On the right side of the editor a list of commands (Steps) of which the program is composed of is shown. After selecting a program command from the list in the left window it can be inserted into the desired place using the  button.
4. Using a double click or the 'Edit' button  parameters, single steps can be changed later.
5. Selected program steps can be deleted from the program by using the 'Delete' button .
6. Use the 'Save' button  to save a program.
7. Save and start programs in one step by simply using the 'Save&Run' button .
8. For even faster programming the VWR thermal cycler comes with a selection of template programs, which can easily be edited and saved under a new program name.

Adding and removing program steps.

Program name.

Program window.

Available commands:

Lock lid (except for Ristretto)
Locks the lid avoiding accidental opening. By activating Lid Pressure On, you can set lid pressure on the sample between 100 - 250 N. When using a HPL Thermal cycler.

Cycle
Combines Open Loop, Close Loop, Temperature, Ramp and Gradient.

Unlock lid (except for Ristretto)
Unlocks the lid if the lid has been locked in a previous program step.

Pause
Interrupts the program temporarily at defined temperature settings.

Heat lid
Turns the heated lid on/off.

Store
Holds the defined temperature, turns off the heated lid automatically.

Temperature
Temperature step for definable time, Temperature and Time Incr./Decr. function.

Start cycle
Begins manual loop.

Gradient
Realises different temperatures at the same time for protocol optimisation (see 'FlexGradient' below for more information).

Close cycle
Ends manual loop.

Sample program

Create a new program with 'New'. Enter a program name and confirm with **OK** to access the programming functions.

01.

Heat lid

1 Heat Lid to 110.0°C

Heated Lid on, temperature 110 °C

A 'Heat Lid' step is inserted in each program by default as a first step. This step can be edited or deleted later. Standard settings can be changed under 'System' => 'Firmware' => 'Settings'.

For 'Store' or 'Temp' steps <30 °C the lid heating is switched off automatically.

To enter another program step select the desired command from the list of the available commands and add it to the program window using the '+' arrow. As a result the following corresponding parameters will be asked.

02

Temp

2 Temp. 95.0°C for 2' 0"

Denaturation at 95 °C for 2 minutes

03 -
07

Cycle

3 Start Cycle, 30x
4 Temp. 95.0°C for 30"
5 Temp. 65.0°C (1.0°C/s) for 30" (~...
6 Temp. 72.0°C for 1' 0" (+0.0°C,+...
7 Close Cycle

Open the 'Cycle' function

→ Enter number of cycles – 30 cycles

→ 3 x 'Temp' step

→ Terminate 'Cycle' function with 'Close Cycle'

04

Temp

4 Temp. 95.0°C for 30"

Denaturation at 95 °C for 30 seconds

05	 Temp	5 Temp. 65.0°C (1.0°C/s) for 30" ~...	Annealing at 65 °C in the 1st cycle for 30 seconds 'Touchdown' - PCR via Temperature-'Decrement'. The annealing temperature of 65 °C in the 1st cycle is reduced by 0,4 °C each cycle. This leads to high specificity in the first cycles and high amplification efficiency in the later cycles. The reduction of the cooling rate at this step of the program is not necessary, but should only be an example. It might make sense, for example, for quantitative annealing of complementary oligo nucleotide or at annealing steps in the first cycles, when primer annealing should be very efficient because of low template amount (e.g. after cDNA synthesis).
06	 Temp	6 Temp. 72.0°C for 1' 0" (+0.0°C,+...	Elongation at 72 °C for 1 minute in the 1st cycle Via Time -'Increment'- the elongation time is prolonged by 5 s each cycle. The incremental prolonging of the elongation time is useful with high numbers of cycles (>30), when polymerase activity declines.
07		7 Close Loop	Terminate 'Cycle' function with 'Close Cycle'.

Alternative programming of steps 03 – 07:


03	 Start Cycle	3 Start Cycle, 30x	Open program cycle – 30 cycles.
04	 Temp	4 Temp. 95.0°C for 30"	Denaturation at 95 °C for 30 seconds
05	 Temp	5 Temp. 65.0°C (1.0°C/s) for 30" ~...	Annealing at 65 °C in the 1st cycle for 30 seconds 'Touchdown' - PCR via Temperature - 'Decrement'. The annealing temperature of 65 °C in the 1st cycle is reduced by 0,4 °C each cycle. This leads to high specificity in the first cycles and high amplification efficiency in the later cycles. The reduction of the cooling rate at this step of the program is not necessary, but should only be an example. It might make sense, for example, for quantitative annealing of complementary oligo nucleotide or at annealing steps in the first cycles, when primer annealing should be very efficient because of the low template amount (e.g. after cDNA synthesis).
06	 Temp	6 Temp. 72.0°C for 1' 0" (+0.0°C,+...	Elongation at 72 °C for 1 minute in the 1st cycle Via Time - 'Increment' the elongation time is prolonged by 5 s each cycle. The incremental prolonging of the elongation time is useful with high numbers of cycles (>30), when polymerase activity declines.
07	 Close Cycle	7 Close Cycle	Terminate program loop with 'Close Cycle'.
08	 Temp	8 Temp. 72.0°C for 5' 0"	Final elongation at 72 °C for 5 minutes The final elongation ensures the complete synthesis of all amplifications.
09	 Store	9 Store forever at 8.0°C	8 °C forever, lid heat is turned off At the 'Store' step the lid heat is turned off automatically. This prevents the lid heat and Peltier elements of the block working against each other. This can happen especially at longer incubation times at low temperature – doing this will increase the life time of the Peltier elements and reduce the power consumption.

FlexGradient Technology (except for Ristretto):


The temperature of each of the eight rows of the block can be individually controlled! For the first time you can choose either the perfectly linear temperature gradient, ideally suited for PCR optimisation, or independent lane control for the use of different primer pairs in one PCR run.

The temperature range is defined by the 'Center' and 'Gradient $\pm$ ' buttons. The standard gradient is always shown from the top to the bottom over the 8 rows of the block. Moreover, with the UNO⁹⁶ you have the opportunity to set 16 gradients. In this case different temperatures range from left to right in the thermo block.

You will also have a gradient in the lanes from left to right - the temperatures displayed in the software are only shown in rows 3 and 10.

By selecting one of the temperature buttons **57.9** different options are available: When the **linear** button is deactivated beforehand, the single temperature of each row can be adjusted within a special range using the 'Adjust Temperature' button. In addition a selected temperature can be used as a new centre value or the gradient step can be changed in a temperature step using the selected temperature.

Program options

When setting up a new program or when saving during programming (via the Info button) program options can be set. The following options are available for the user:

1. Power Fail Denaturation

If a power failure happens while a program is running the thermal cycler will restart the program at the point it was stopped as soon as the power returns. If the 'Power Fail Denaturation' option is activated, the thermal block will heat to the stated temperature for the set time before continuing with the program.

2. Global Program Ramp

By using this function a global ramp rate in °C/s can be programmed. This global ramp rate will be used in all temperature and gradient steps of the PCR program, but can be changed in a single temperature or gradient step manually by entering the desired ramp value at the step.

3. Gradient Control

Using this option you can decide whether in a gradient step the plateaus (plateaus synchronised) or the ramps (ramps synchronised), should be synchronised for all rows.

4. Tube Control


Using tube control gives you the ability to adjust temperature regulation for a temperature step according to sample volume. If necessary, the temperature regulation will perform short-term deliberate overshoots and undershoots depending on the entered sample volume during PCR run. Especially if combined with the emulation mode, this function can be helpful if the emulated cycler uses the tube control function before too.

5. Emulation


Using the emulation function the thermal cycler will perform PCR similar to a chosen cycler, emulating its performance. This function could be useful when transferring a PCR program to the UNO or Doppio from another type and model of thermal cycler.

5.5 Further dialogues: Diagnostics, GLPs, System


Diagnostics: Monitoring the actual temperatures of each of the six Peltier elements during a run. In the event of an instrument malfunction have a look at the Peltiers' temperature profile.

GLPs: Displaying and organising GLP (Good Laboratory Practice) reports created by the thermal cycler. This option is mainly for labs that need to save GLP reports for documentation purposes.

System: Configuration of the system and the general changing device settings.

Function buttons

	Time Setting of date and time.		Hide/show cursor Hide or shows the cursor on the screen.
	User call Setting of acoustic signals or email notification (Pause, Store, Program End).		User management Defines user and user groups (plus their access rights). The default password for the user 'Administrator' is: 'admin' . Via 'Logout' button the current user can be logged off. As a result the user 'Guest' is automatically logged in. This user will only have minimal user rights. These settings can be changed by the administrator if desired.
	Network settings For network integration of the VWR Cyclar.		
	Hardware Displays hardware versions of Control Unit, Power Unit and Thermo block.		Firmware Displays, administrates and updates software (Control/Power Unit, Thermo block).
	Printer Sets up a printer path for an existing network printer.		Display settings Changes settings of the screen display.
	Protocol Internal communication of the thermal cyclar for failure diagnosis.		Service Service function (only for authorised technical service)
	LAN control For PC and Master/Slave control.		USB recovery Manual mounting of USB sticks that are not automatically detected.

5.6 Further questions...?

Direct help :

Use this button to activate (deactivate) the direct help mode. If the mouse pointer is visible a question mark appears alongside. By touching a function button or other dialogue element the corresponding help window will appear.

A more detailed instruction manual in English or German is available for download on our website vwr.com.

Technical service and ordering information

Should you have any questions about the VWR thermocycler, please visit vwr.com to find the respective contact person.

This short manual is delivered with the following products:

Description	Delivered with	Cat. No.
Doppio	2x 48-well Universal blocks with standard lids for 48x 0,2 ml tubes, 48-well PCR plates or 24x 0,5 ml tubes with flat caps each	732-2551
Doppio Gradient	2x 48-well Universal Gradient blocks and standard lids for 48x 0,2 ml tubes, 48-well PCR plates or 24x 0,5 ml tubes with flat caps each	732-2552
UNO ⁹⁶	96-well Universal block and standard lid for 96x 0,2 ml tubes, 96-well PCR plates or 48x 0,5 ml tubes with flat caps	732-2548
UNO 96 HPL	96 well Universal block with High Pressure lid (HPL, 100 – 250 N) for the safe sealing of 96 x 0.2 ml tubes, 96 well PCR plates or 48 x 0.5 ml tubes with flat caps.	732-2914
UNO ⁹⁶ Gradient	96-well Universal Gradient block and standard lid for 96x 0,2 ml tubes, 96-well PCR plates or 48x 0,5 ml tubes with flat caps	732-2549
UNO 96Gradient HPL	96 well Universal block with High Pressure lid (HPL, 100 – 250 N) for the safe sealing of 96 x 0.2 ml tubes, 96 well PCR plates or 48 x 0.5 ml tubes with flat caps	732-2915
UNO ³⁸⁴	384 well block with High Pressure Lid	732-2550
Ristretto	32 well Universal block with lid heating and automatic height adaption for 32 x 0.2 ml tubes or 16 x 0.5 ml tubes with flat caps	732-2553

6 TECHNICAL DATA

General characteristics

- 8, 2 x 8 or 16 Peltier elements with Long-Life-technology and 8, 2 x 8 or 16 control circuits with Pt 1000 temperature sensors, respectively (at Ristretto, Doppio and UNO)
- Maximum heating and cooling rate: 5 °C/s (3 °C/s for Ristretto)
- Block uniformity (at 72 °C): $\pm 0,20$ °C
- Temperature range thermo block: 4 to 105 °C
- Control accuracy thermo block: $\pm 0,1$ °C
- Variable ramping: 0,1 to 3,0 °C/s
- Increment/decrement time: 0:01 to 9:59 min
- Increment/decrement temperature: 0,1 to 9,9 °C
- Lid heating temperature range: 40 to 120 °C
- Programmable lid locking mechanism
- 4 x USB (1x USB for Ristretto), 1x Ethernet interfaces (MS Windows® or Linux), and maximum length of power cable is 3 m!

Note: USB interfaces only support standard USB sticks, cable mouse and keyboards!

- Dimension (WxHxD): 30x28x38 cm, Ristretto: 27x19x33 cm
- Weight including block: 13,3 kg (5.2 kg for Ristretto)
- Electrical power supply: 100 – 240 V AC, 50/60 Hz, 850 VA
- Fuse: 10 A time lag
- Degree of pollution: 2
- Ambient temperature: 10 to 30 °C
- Maximum relative humidity: 70%
- Maximum altitude for the use of this system: Up to 2000 m above sea level
- For indoor use only!
- Noise level: <37 dB (A)

Gradient features (optional)

- Maximum gradient range over 8 or 16 columns, respectively: 30 °C (± 15 °C)
- Temperature range gradient: 35 to 105 °C
- Gradient accuracy: $\pm 0,1$ °C
- Gradient function can be upgraded even after purchase and on-site
- FlexGradient Technology: Ideal, linear gradient or independent lane control
- 'Gradient Control': Synchronised plateau times or ramp rates

User interface/functions

- Touch-sensitive TFT-Display (VGA, Graphic, 65535 colours), operation via USB mouse possible
- Direct help function, T_m- and elongation time calculator
- Internal flash memory for 500 000 typical PCR programs in free configurable folders/subfolders
- Unlimited number of programs via network PC or USB memory stick
- Free computer software for the generation of PCR programs on the computer
- Unlimited number of steps/programs
- Password protected user accounts with variable access rights
- GLP reports for complete recording of all runs
- Quick start function of the last used program
- Automatic restart after power failure, Power-Fail-Denaturation, Instant Incubation
- Remote control of the system by PC software
- Free choice of mp3 signal sounds
- User-notification by email
- Master/Slave control (wireless option)

Block options (partly available as interchangeable blocks)

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2x 48-well Universal blocks with lid heating and automatic height adaption for 48x 0,2 ml tubes, 48-well PCR plates or 24x 0,5 ml tubes with flat caps.

VWR UNO⁹⁶/96 Gradient:

96-well Universal block with lid heating and automatic height adaption for 96x 0,2 ml tubes, 96-well PCR plates or 48x 0,5 ml tubes with flat caps.

UNO 96 HPL /96 HPL Gradient:

96-well Universal block with High Pressure lid (HPL, 100 – 250 N) for the safe sealing of 96 x 0.2 ml tubes, 96-well PCR plates or 48 x 0.5 ml tubes with flat caps.

VWR UNO³⁸⁴:

384-well block with High Pressure Lid (HPL, 100 – 250 N) for the safe sealing of 384-well PCR plates.

Ristretto:

32 well Universal block with lid heating and automatic height adaption for 32 x 0.2 ml tubes or 16 x 0.5 ml tubes with flat caps

7 MAINTENANCE AND REPAIR OF THE THERMAL CYCLER

7.1 Cleaning, maintenance and repair

For safety reasons the device must be switched off and the mains plug pulled out before cleaning is carried out! The surface of the thermal cycler can be wiped off with a damp cloth or ethanol. Do not use aggressive or scouring cleaners or organic solvents for cleaning. The device should be protected from aggressive chemicals. Make sure that no liquid penetrates the interior of the device.

If sample liquids penetrate into the drill holes of the thermal cycler clean them immediately with mild soap followed by distilled water. You can use mild disinfectants for decontamination.

Maintenance (over temperature test and temperature check up, see 8.2 and 8.3) should be carried out regularly to detect possible technical failures.

The electric fuses are situated at the back of the device between the power switch and the mains plug. After pulling out the mains plug and pushing the safety catch upwards the fuses can be changed. Use only fuses with correct values (indications about the fuse type are situated at the back of the device).

The device should only be opened by a qualified specialist. Unauthorised work on the device voids the warranty. All repairs should only be carried out by authorised persons using original replacement parts. The replacement of single components of the device should not be done by the user but exclusively by authorised specialists. Please therefore send the device back to the manufacturer. The surface of the device should be decontaminated with a lint-free cloth soaked with 70% ethanol before sending. The drill holes of the thermo block should be cleaned with a 1,5% hypochlorite solution with the help of a cotton bud. Please fill in and sign the decontamination certificate and send it to the manufacturer together with the device.

7.2 Over temperature test

A periodical test of the over temperature protection (→ System → Hardware → Thermo Block → Over temperature test) is recommended.

7.3 Temperature check up

Maintenance of the temperature at regular intervals (e.g. annually) is recommended. The thermometer should be equipped with wired sensing devices (e.g. PT1000) and should have an accuracy of 1/10 class B ($\Delta T = \pm 0,03 \text{ }^{\circ}\text{C}$) measurements in the range of 0 to $+100 \text{ }^{\circ}\text{C}$. While processing the temperature control room temperature should not exceed $+25 \text{ }^{\circ}\text{C}$.

7.4 Calibration

If there is a reasonable suspicion that the device has deviated from its calibration range, because the variations of the temperature check up are considerably out of the acceptable tolerance, a new calibration must be carried out by the manufacturer. Please contact your local distributor.

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3 PRESENTATION DU SYSTEME

VWR thermal cyclers are the ideal PCR devices for all molecular biological and biochemical laboratories working in the field of basic research or routine diagnostic. The 250 °C HTR (High Temperature Range) Peltier technology allows quick changes in temperature up to 5 °C/sec (max). The individual monitoring and controlling of each of the eight or 16 Peltier elements per block equipped with 'Long Life Technology' ensure the tightest correlation of target and actual temperature across the entire thermal plate leading to highest reproducibility of the results. To avoid condensation each block is equipped with a heated lid. The VWR thermal cyclers are equipped with a touch-sensitive, graphic coloured TFT display of the newest generation, allowing operation of the device by simply touching symbols and graphic elements on the screen. This leads to easy and intuitive handling and programming of the device. Additionally a mouse can be connected via USB (reboot required). Furthermore an Ethernet port is available to connect the devices with Microsoft Windows® or Linux networks. It is therefore possible to use an external server for archiving programs and GLP reports which will be available globally in the network (e.g. for other thermal cyclers). Transfer of data can be done by using commercially available USB memory sticks enabling a nearly unlimited saving space for programs and GLP reports, which can be printed using a printer connected to the thermal cycler via the network.

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

The screenshot shows a 'Gradient' control window. It has two main sections: '8-Gradient' and '16-Gradient'. The '16-Gradient' section is active, showing a linear gradient from 55.0°C to 65.0°C. The 'Center' temperature is set to 60.0°C, and the 'Gradient ±' is 5.0°C. The 'Time' is 0 h 0 min 0 s, and the 'Ramp Rate' is 0.0°C. The 'Temp.-Incr./Decr.' is 0.0°C, and the 'Time-Incr./Decr.' is 0 min 0 s. The interface includes a 'linear' checkbox and a 'convert this Gradient Step to a Temperature Step with 57.9°C?' button. The bottom of the window shows 'Row 3' and 'Row 10' with corresponding temperature values.

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The screenshot shows a dialog box with three buttons: 'Adjust Temperature', 'use 57.9 as a new Center Temperature', and 'convert this Gradient Step to a Temperature Step with 57.9°C?'. A red 'X' button is also visible.

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VWR International AB
VWR International AG
Teknolojileri Ltd.Şti.
VWR International Ltd