
Digitrate *Pro*

Instruction Manual

For capacities of 30 ml and 50 ml

	EU Cat. No.	NA Cat. No.
30 ml	612-3875	75856-476
50 ml	613-5287	75856-478
Version: 1		
Issued: May 2017		

DIGITRATE *PRO*

 EN	Instruction Manual.....	1 – 14
 FR	Manuel d'instructions	15 – 28
 DE	Bedienungsanleitung	29 – 42
 ES	Manual de instrucciones	43 – 56
 IT	Manuale d'istruzioni	57 – 70
 PT	Manual de instruções	71 – 84

United States

Manufactured for:

VWR International, LLC
100 Matsonford Rd
Radnor, PA 19087
800-932-5000

<http://www.vwr.com>

Europe

Address of legal manufacturer

VWR International bvba
Researchpark Haasrode 2020
Geldenaaksebaan 464
B-3001 Leuven
+ 32 16 385011

<http://be.vwr.com>

Contents

Warning	3
Safety Information.....	4
Package Contents.....	4
Installation.....	5
Restrictions of Use.....	5
Before Using the Digitrate Pro	5
Assembly.....	5
Specification.....	5
Instructions for Use	5
Priming.....	6
Titrating.....	6
Freeze Count Feature.....	6
Auto Off.....	6
Adjustable Head	6
User Calibration Procedure	7
Troubleshooting	8
General Maintenance	9
Maintenance/Cleaning.....	9
Sterilisation/Autoclaving	9
Fluid-Path Disassembly Procedure.....	9
Further Disassembly Options	10
Accessories and Spares.....	12
Adaptors	12
Reservoirs	12
Spare Parts	13
Technical Service	13
Return for Repair	13
Servicing.....	13
Warranty	13
Equipment Disposal.....	14

WARNING!

- Always observe general safety regulations (e.g. wear suitable eye, face, hand and general protection).
- Follow this operating instruction manual.
- NEVER use force! Breakage of any part may lead to the operator and/or other persons being exposed to hazardous substances.
- Do not use the instrument where flammable fumes may exist.
- NEVER use damaged or deformed components.

DIGITRATE *PRO*

Safety Information

Careful safety precautions must apply when dispensing corrosive, radioactive or hazardous chemicals.

- When using the instrument, take into consideration the chemical compatibility of the materials in contact with the liquid.
- Always check the Digitrate *Pro* for leaks. Make sure that the inlet tube (14) and FEP dispense tube (12) are firmly in position before starting to use the instrument.
- Never carry the instrument while holding onto its dispensing barrel or electronic head (2). Always support the bottle with one hand under its reservoir and the other hand around the threaded platform base (10) of the dispenser in order to stabilize the unit.
- Avoid contact with the FEP dispense tube nozzle (12).
- When dispensing make sure that anti-drip valve is open (11).
- Use only original spare parts and original accessories (See Page 12 and 13).
- If in doubt, consult your safety officer.

Package Contents

- Digitrate *Pro* titrator
- Three adaptors 33 mm, 38 mm and 45 mm (See Page 12).
- PTFE inlet tube (See Pages 11 and 13).

Note

Refer to pages 10 and 11 for diagram references.

Installation

Restrictions of Use

NEVER use the Digitrate *Pro* with:

- Liquids which are not compatible with PTFE, PVDF, borosilicate glass, alumina ceramic, FEP or those which might attack platinum-iridium
- Hydrofluoric acid
- Liquids which contain solid particles
- Temperature limits of Digitrate *Pro* and reagent are 15 °C to 40 °C

Before Using the Digitrate *Pro*

Check that the instrument has not been damaged in transit. Do not grip the electronic head casing (2) when screwing, unscrewing or adjusting the instrument.

The Digitrate *Pro* is fitted with a non-rechargeable lithium battery. The battery life is greater than 60,000 three minute titrations (low battery is indicated with an icon on the display). The battery should be changed only by a certified person. The lithium battery should be disposed of correctly, on no account should it be incinerated.

Assembly

The Digitrate *Pro* is packed with the FEP dispense tube (12) attached and the PTFE inlet feed tube removed (14). The 300 mm length of PTFE inlet tube provided should be trimmed to fit your particular container. The bottom end should be cut diagonally at a length that is close to the bottom of the container when the Digitrate *Pro* is assembled onto the container. Longer lengths of inlet tube are available on request (14).

Four adaptors are available as accessories to suit containers with a 33 mm, 38 mm, 40 mm or 45 mm screw neck. (3 are included with this unit). The threaded platform base on the Digitrate *Pro* has a 30 mm screw thread for use with containers of this size.

The assembled Digitrate *Pro* is screwed to the reservoir using gentle hand torque applied to the threaded platform base (10) only. Removal should also be by means of hand torque applied to the threaded platform base (10) only.

Do not operate the titrate wheels (1) until the unit is safely and fully mounted onto the reservoir.

Specification

Accuracy $\leq \pm 0.2\%$ on maximum delivery and a precision $\leq 0.1\%CV$ using distilled water at 20 °C.

Testing takes place according to DIN EN ISO 8655 using maximum nominal volume of distilled water at 20 °C, with a uniform and smooth dispense action.

Instructions for Use

Four momentary action push buttons

1. 'On' - switches the unit on and displays the last reading (15).
2. 'Reset' - reverts display to 00.00 (18).
3. 'CAL' - user custom calibration - pen point access (17).
4. 'Mode' - toggles between aspirate or titrate on the LCD display (16).

DIGITRATE *PRO*

Priming

Position an empty container under the right angle of the FEP dispense tube nozzle (12), open anti-drip valve (11) and move the piston to the fully down position. The FEP dispense tube nozzle (12) must point away from the user at all times. Prime the unit by a few gentle up and down piston strokes, taking the piston right down and lifting about an inch or so. Repeat until a steady, bubble free flow is obtained from the FEP dispense tube nozzle (12). Priming can be completed without switching the unit on. After priming, the piston should be left in the fully down position. When not in use turn anti-drip valve (11) to the closed position.

Titrating

To titrate, first switch the instrument on by pressing the 'ON' button (15). Press 'MODE' (16) until the LCD arrow points upwards. Press 'RESET' (18) to zero the unit. Open anti-drip valve (11) and gently turn the titrate wheels (1) 'up' so filling the borosilicate glass barrel (8). The amount of fluid taken in should be greater than the volume required for titration and will be displayed on the LCD screen.

Note: Too fast a filling action may lead to air bubbles forming under the piston.

After each filling, compensate for backlash by turning the hand wheels (1) slightly in the discharge direction until 1-2 drops are dispensed into an empty container.

To start titration, press the 'MODE' button (16) so that the arrow points down. Press the 'RESET' button (18) to zero the display. Place a receiving vessel under the FEP dispense tube nozzle (12) and gently turn the titrate wheels (1) 'down' so dispensing liquid. The dispensed volume will be shown on the display. Do not drive the piston to the bottom of its stroke, but stop a few ml short of it. If necessary, the borosilicate glass barrel (8) may be refilled at any time during the titration, without affecting the dispensed volume shown on the LCD using 'Freeze Count' feature.

Freeze Count Feature

Digitrate *Pro* has a unique 'Freeze Count' feature enabling the user to freeze the displayed value. This function allows backlash compensation after replenishing the barrel during titration, without affecting the value displayed. Simply press and hold 'ON' button (15), aspirate to refill barrel. Release 'ON' button (15) to continue titration. It also allows the user to take a sample of the dispensed liquid during titration, just press 'ON' button (15) to 'Freeze Count'.

Auto Off

The Digitrate *Pro* will automatically switch off when it is idle for a period of more than three minutes. The instrument is switched on by pressing the 'ON' switch (15). When the instrument is switched on, the previous value is displayed. For subsequent titration, press the 'RESET' button (18) to zero the display.

Adjustable Head

The display unit head (2) can rotate on the threaded platform base and be adjusted to suit the user.

Adjust the orientation of the display head (2) to the desired position, only when the Digitrate *Pro* barrel is empty. Loosen the threaded ring (4) of the platform base by a quarter turn and rotate the display head (2) to the desired angle by holding the threaded base (10) of the Digitrate *Pro* with one hand, while the other hand rotates the polypropylene protection sleeve (3). When the angle of the display is in its desired position, retighten the threaded ring.

Never rotate the display by holding on to the 'electronic display head casing' (2).

User Calibration Procedure

The user is allowed to override Factory Calibration and custom calibrate the Digitrate *Pro*.

The procedure is as follows:

1. Establish a mean volume of 10 dispensings of nominal volume of distilled H₂O at 20 °C (e.g. 50 ml for the 50 ml unit), and establish the mean value gravimetrically. After each filling, remember to reset the display (18) to 00.00, only after you have compensated for any backlash.
2. When the mean volume, taking into account the correction factor (see formula below), has been established, turn the titrate wheels so that nominal volume is displayed i.e. 50 ml for the 50 ml unit.
3. With the instrument turned on, press the calibration switch (17) through the right hand pen hole.
4. The LCD will flash 'CAL'.
5. Use the 'ON' button (15) as '-' and the 'Reset' button (18) as '+' to set the calculated mean volume on the display. Push once for increment of '1' – continuous push runs the numbers faster.
6. Store the new value by pressing the 'MODE' button (16). The word 'CAL' appears on the display from now on in order to indicate that the unit has been calibrated by the user.

Calculation (for the nominal volume)

$$\text{Mean value } \bar{x} = \frac{\sum x_i}{n}$$

x_i = results of the weighings

n = number of weighings

$$\text{Mean volume } V = \bar{x} \cdot Z$$

Z = correction factor (e.g. 1.0029 µl /mg at 20 °C, 1013 hPa)

Calculate the dispensed volume by taking the temperature into account
(1 ml distilled water at 20 °C = 0.997 g)

Factory calibration figures, however, are retained in the memory at all times.

To return to the factory calibration:

1. Press the calibration button through the pen hole twice.
2. The display cancels the word 'CAL' and the unit is returned to the factory calibration.

Troubleshooting

Problem	Possible Cause	Remedy
Air bubbles appear in discharge nozzle	Liquid reservoir is empty	Refill reservoir and prime unit
	Too fast filling action	Operate titrate wheels (1) smoothly/slowly to ensure all air is removed
	Glass barrel (8) is not sealing against FEP 'O' ring (9)	Unscrew the threaded platform ring (4) and make sure that the FEP 'O' ring (9) is properly fitted into its recess
	Leaking piston (7)	Clean PTFE piston (7). If problem persists, replace piston
	Leaking discharge valve	Clean unit by flushing through - if problem persists, replace platform base (10)
Barrel does not fill with liquid	Intake tube not fitted correctly (14)	Connect correctly
	Inlet valve stuck (11)	Free inlet valve (11) by inserting a thin rod into the inlet aperture and push gently
	Glass barrel (8) is not sealing against the FEP 'O' ring (9)	Unscrew the threaded platform ring (4) and make sure that the FEP 'O' ring (9) is properly fitted into its recess
Titrating not possible	Blocked dispense nozzle (12)	Disassemble the dispense nozzle (12) and flush through with cleaning fluid
	Discharge valve stuck (11)	Clean unit by immersing platform in cleaning fluid - if problem persists, replace platform base (10)
Wrong dispense volume	Instrument not calibrated	Calibrate the unit
	Backlash not compensated for	Compensate for backlash
	Leaking valves (11)	Clean platform base (10) - if problem persists, replace
	Titrate wheels (1) not turning uniformly	Operate titrate wheels (1) slowly
Liquid appears between glass barrel and polypropylene protection sleeve	Glass barrel (8) is not sealing against the FEP 'O' ring (9)	Unscrew the threaded platform (4) ring and make sure that the FEP 'O' ring (9) is properly fitted into its recess, slightly stretching the 'O' ring with your fingertips may assist this
	Damaged FEP 'O' ring (9)	Replace FEP 'O' ring (9)
Blank display despite 'ON' status	Empty battery	Send unit for battery replacement
	Malfunctioning electronics	Send unit for repair

General Maintenance

Maintenance/Cleaning

1. Place the instrument into a sink together with its reservoir. Unscrew the threaded platform base (10) and lift the dispenser's intake tube (14) carefully out of the reservoir, whilst tapping it against the reservoir's aperture in order to shake off any droplets from the intake tube.
2. Hold the FEP dispense tube nozzle (12) over the aperture of the reservoir and apply gentle piston strokes in order to return the contents into the reservoir.
3. Flushing out with distilled water or a suitable solvent is recommended after the Digitrate *Pro* has been used with crystallizing fluids. This will preserve the smooth action of the piston and free action of the inlet and outlet valves. If the inlet valve (11) does stick and is not freed by flushing, it may be freed by gently inserting a thin rod into the inlet aperture and gently pushing the ball off its seating. Re-check unit for operation.
4. Empty the instrument completely after cleaning.

Note:

All maintenance should be carried out wearing suitable eye protection and protective clothing. If in doubt, consult your safety officer.

Sterilisation/Autoclaving

Note:

The electronic display head cannot be sterilised or autoclaved. Chemical sterilisation will damage electronic components. Autoclaving will cause permanent damage to the electronic components and the battery may explode.

Before sterilisation or autoclaving, please follow the maintenance 'flushing out procedure'. Sterilisation/autoclaving is carried out to 'fluid-path' components only.

Chemical sterilisation is permissible by soaking the 'fluid-path' components overnight in a dilute (1:1000) solution of sodium hypochlorite. The Digitrate *Pro* should be rinsed well in sterile, distilled water if the latter method is used.

Autoclaving of 'fluid-path' components is permissible at 121 °C, 2 bar, after the normal cleaning procedure has been carried out.

Fluid-Path Disassembly Procedure

Disassembly should only be undertaken AFTER the unit has been cleaned, using the standard cleaning procedure described above. Undo the threaded platform ring (4) of the base platform and remove the electronic display head (2) together with the PTFE piston (7), piston rack (5) and polypropylene protection sleeve (3) from the platform base (10). Take care not to lose the FEP 'O' ring (9) which is mounted inside the platform recess.

Slide the glass barrel (8) off the piston (7) and lay it in a safe place. Turn the titrate wheels (1) down until the piston rack (5), together with the PTFE piston (7), disengages from the electronic display head (2) and slide them back into the glass barrel (8) from the bevelled edge end, taking care not to damage the piston when inserting it into the glass barrel.

Place these components and the threaded base platform (10) with FEP dispense tube (12) and protection cover (13) attached, together with PTFE inlet tube (14) on a suitably soft surface within the autoclave in order to avoid metal contact.

Steam sterilisation of the piston outside the glass barrel may damage it. Do not use hard tools to scrape off residue of reagent from the piston. Do not use force in assembly or disassembly.

DIGITRATE *PRO*

Further Disassembly Options

To clean or replace the FEP dispense tube (12) and/or threaded platform base (10), follow these steps;

Turn the anti-drip valve (11) to the closed position and unscrew the threaded ring (4) to separate the pedestal base from the glass barrel of the Digitrate *Pro*. Remove dispense tube protection cover (13) as shown opposite in an upward direction allowing FEP dispense tube (12) to be pulled out of anti-drip valve assembly (11) in direction indicated.



Reassemble in reverse, make sure the front end of the FEP dispense tube (12) is clipped into the dispense tube protection cover (13) from the underside.



Reassemble the remaining components in reverse order, making sure that the FEP 'O' ring (9) is well located into its platform recess.

Note:

Bevelled inside edge of barrel must be at the top end when fitted. The piston rack must engage in to the electronic head with the piston rack teeth pointing forward in the direction of the display.

After reassembling the instrument, prime with distilled water to ensure that assembly has been correctly followed and piston is working smoothly. Check that no leaks occur.


Accessories and Spares


Adaptors

A range of four adaptors for fitting your dispenser to the reservoir are available along with an angled funnel to aid filling of reservoirs.


33 mm

EU Cat. No. 612-3980
NA Cat. No. 53409-990


38 mm

EU Cat. No. 612-3978
NA Cat. No. NA 53409-992

40 mm


EU Cat. No. 612-3977
NA Cat. No. 53409-994

30 mm

45 mm

EU Cat. No. 612-3979
NA Cat. No. 53409-996

115 mm

Angled
Funnel

EU Cat. No. 613-0969
NA Cat. No. 14526-004

Reservoirs

Description	EU Cat. No.	NA Cat. No.
1.4 litre 2-neck, amber Squat	613-0971	14526-000
2.5 litre 2-neck, amber Squat	631-0975	14526-002
2.5 litre 2-neck, amber Winchester	613-0972	53410-044

Spare Parts

Description	Drawing	30 ml		50 ml	
		EU Cat. No.	NA Cat. No.	EU Cat. No.	NA Cat. No.
Titrate Wheels	1	613-5653	10782-760	613-5653	10782-760
Threaded Platform Ring	4	613-5663	10782-852	613-5664	10782-854
Borosilicate Glass Barrel	8	613-5655	10782-766	613-5654	10782-764
FEP 'O' Ring	9	613-5652	10782-758	613-5651	10782-756
FEP Dispense Tube	12	613-5659	10782-790	613-5659	10782-790
Polypropylene Dispense Tube Protective Cover	13	613-5660	10782-794	613-5660	10782-794
Inlet Tube	14	613-5658	10782-772	613-5658	10782-772
Piston & Rack Complete	19	613-5657	10782-770	613-5656	10782-768
Digitrate <i>Pro</i> Pedestal Assembly	20	613-5661	10782-812	613-5662	35001-360

See page 10/11 for reference

Technical Service

Return for Repair

In the unlikely event of repair, or where the damage to the unit necessitates return to your distributor or the manufacturer, please decontaminate the unit and complete the required decontamination certificate where appropriate. In addition, please give written details of the fault. The decontamination certificate is available from the product ordering page on www.vwr.com.

Servicing

The manufacturer advises servicing this unit annually to prolong the life of the unit. For more information please contact your local VWR representative or go to the product ordering page on www.vwr.com.

WARRANTY

The manufacturer warrants that this product will be free from defects in material and workmanship for a period of three (3) years from date of delivery. If a defect is present, VWR will, at its option and cost, repair, replace, or refund the purchase price of this product to the customer, provided it is returned during the warranty period. This warranty does not apply if the product has been damaged by accident, abuse, misuse, or misapplication, or from ordinary wear and tear. If the required maintenance and inspection services are not performed according to the manuals and any local regulations, such warranty turns invalid, except to the extent, the defect of the product is not due to such non performance.

Items being returned must be insured by the customer against possible damage or loss. This warranty shall be limited to the aforementioned remedies. IT IS EXPRESSLY AGREED THAT THIS WARRANTY WILL BE IN LIEU OF ALL WARRANTIES OF FITNESS AND IN LIEU OF THE WARRANTY OF MERCHANTABILITY.


Equipment disposal

This equipment is marked with the crossed out wheeled bin symbol to indicate that this equipment must not be disposed of with unsorted waste.

Instead it is your responsibility to correctly dispose of your equipment at lifecycle-end by handing it over to an authorized facility for separate collection and recycling. It is also your responsibility to decontaminate the equipment in case of biological, chemical and/or radiological contamination, so as to protect from health hazards the persons involved in the disposal and recycling of the equipment.

For more information about where you can drop off your waste equipment, please contact your local dealer from whom you originally purchased this equipment.

By doing so, you will help to conserve natural and environmental resources and you will ensure that your equipment is recycled in a manner that protects human health.

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
GO TO **VWR.COM** FOR THE LATEST NEWS,
SPECIAL OFFERS AND DETAILS OF YOUR
LOCAL VWR DISTRIBUTOR

Prices and product details are current when published; subject to change without notice. | Certain products may be limited by country, federal, state, provincial, or local regulations. | VWR makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC. Offers valid in countries listed above, void where prohibited by law or company policy, while supplies last. | VWR, the VWR logo and variations on the foregoing are registered (®) or unregistered trademarks and service marks, of VWR International, LLC and its related companies.

All other marks referenced are registered by their respective owner(s). | Visit vwr.com to view our privacy policy, trademark owners and additional disclaimers.

©2017 VWR International, LLC.

All rights reserved.