

# EVE™ PLUS

The World's Fastest Automated Cell Counter

# LESS THAN 1 SEC!




**FAST**



**EASY**



**ACCURATE**



System Platform	Tablet PC	Data Export	PDF , CSV	Cell size range	5 ~ 60 um
Time of analysis	< 1 sec (Manual Focus) < 10 sec (Auto Focus)	Connectivity	WiFi	Loading volume	10 uL
		Measurement range	1x10 <sup>4</sup> ~ 2x10 <sup>7</sup>	Staining method	Trypan blue stain

#### NanoEnTek, Inc.


851-14, Seohae-ro, Paltan-myeon, Hwaseong-si, Gyeonggi-do,  
18531, Korea Tel : +82-2-6220-7940 / Fax : +82-2-6220-7999

#### NanoEnTek America, Inc.


240 Bear Hill Road, Suite 101, Waltham, MA 02451, USA  
Tel : +1-781-472-2558 / Fax : +1-781-790-5649

[website](http://www.nanoentek.com) www.nanoentek.com  
[e-mail](mailto:sales@nanoentek.com) sales@nanoentek.com

## Correlation of EVE™ PLUS and manual counting





Measuring from the EVE™ PLUS extends further the high concentration range than hemocytometer readings



EVE™ PLUS measures the number and concentration of total cells, viable cells, dead cells. It provides viability results and cell size gating.

## Total cell counting result (Clumped cell)



Clumped cells were counted with EVE™ PLUS, nucleus counting, and the automated cell counters A, B, C. NanoEnTek's ADAM-MC2 device was used for nucleus counting. It is accurate in clumped cells by counting the stained nuclei through the PI staining method. EVE™ PLUS is comparable to the nucleus counting for all cell lines with accuracy and precision. The other automated cell counters A, B, and C were shown inaccurate numbers in the clumped cell. EVE™ PLUS identifies and counts the individual cells within the clumpy cells for accurate analysis.

## Cell lines validated on EVE™ PLUS

Cell Type	Animal	Organ	Growth Properties
HeLa	Human	Skin	Adherent
NIH-3T3	Mouse	Embryo	Adherent
U-2 OS	Human	Bone	Adherent
Jurkat	Human	Blood	Suspension
KG-1	Human	Blood	Suspension
HepG2	Human	Liver	Adherent
Hep3B	Human	Liver	Adherent
LNCaP	Human	Prostate	Adherent
SH-SY5Y	Human	Brain	Adherent
SCN2.2	Rat	Brain	Adherent
F9	Mouse	Embryo	Adherent
MCF7	Human	Breast	Adherent
A549	Human	Lung	Adherent
GH3	Rat	Pituitary gland	Adherent