



Product Specification Sheet

VWR BOTTLES, AMBER, WIDE MOUTH, WITH SCREW CAP

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

Cat. No.	Capacity	Ø ext.	Height	Neck I-Ø	Neck Ø ext.	Pk
215-5676	50 ml	39,0 mm	77,0 mm	24,0 mm	32,0 mm	10
215-5677	100 ml	47,0 mm	94,0 mm	24,0 mm	32,0 mm	10
215-5678	250 ml	63,0 mm	119,0 mm	31,0 mm	40,0 mm	10
215-5679	500 ml	78,0 mm	152,0 mm	39,0 mm	50,0 mm	12
215-5680	1000 ml	96,0 mm	206,0 mm	50,0 mm	65,0 mm	12

Material: Low Density Polyethylene (LDPE)

Other Properties:

- Coloured amber to protect light sensitive substances.
- Leakproof, tough and versatile.
- Perfect for sample collection, storage and transport.
- Ideal for liquids, pastes and granular material.
- Wide mouth allows easy filling and cleaning.

Disclaimer: VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use.

Date of issue: **July 2012**

