

CFC selective supplement

Code 84746.0001

Intended use

Freeze-dried selective supplement for the preparation of CFC medium used for the isolation and enumeration of *Pseudomonas* spp. in meat products (ISO 13720).

Vial contents

Cetrimide..... 5.0 mg
Fucidin..... 5.0 mg
Cephalosporin..... 25.0 mg

for preparation

Prepare 500 ml of *Pseudomonas* agar base (Art. N°: 84650.0500), autoclaved and cooled to 45-50°C. Aseptically reconstitute the contents of one vial of CFC selective supplement with 2 ml of ethanol/sterile purified water (1:1). Mix well and aseptically add to 500 ml cooled medium base. Pour into sterile Petri dishes under aseptic conditions.

Instructions for use

For laboratory use only.

Consult the technical sheet of *Pseudomonas* agar base, Art. N°: 84650.0500

Consult the material safety data sheet before the use.

Physical characteristics

Appearance of freeze-dried supplement	High, white, homogeneous pastille
Appearance of solution after reconstitution	Whitish, turbid
pH of complete medium (25°C)	7.2 ± 0.2

Storage conditions

Store at 2-8°C. When stored as directed the supplement remains stable until the expiry date shown on the label. Do not use beyond stated expiry date.

Ordering information

84746.0001 CFC selective supplement 10 vials, each for 500 ml of complete medium