

PP selective supplement

Code 84747.0001

Intended use

Freeze-dried selective supplement for the preparation of PP medium used for the isolation and enumeration of *Pseudomonas* spp. in milk and milk products (ISO/TS 11059)

Vial contents

Penicillin G.....	50000 IU
Pimaricin.....	5.0 mg

Instructions for preparation

Prepare 500 ml of *Pseudomonas* agar base (Art. N°: 84650.0500), autoclaved and cooled to 45-50°C. Aseptically reconstitute the contents of one vial of PP selective supplement with 2 ml of sterile purified water. Mix well and aseptically add to 500 ml cooled medium base. Pour into sterile Petri dishes under aseptic conditions.

Instructions for use

For laboratory use only.

Consult the technical sheet of *Pseudomonas* agar base, Art. N°: 84650.0500

Consult the material safety data sheet before the use.

Physical characteristics

Appearance of freeze-dried supplement	High, white, homogeneous pastille
Appearance of solution after reconstitution	Whitish, turbid
pH of complete medium (25°C)	7.2 ± 0.2

Storage conditions

Store at 2-8°C. When stored as directed the supplement remains stable until the expiry date shown on the label. Do not use beyond stated expiry date.

Ordering information

84747.0001	PP selective supplement	10 vials, each for 500 ml of complete medium
------------	-------------------------	----------------------------------------------