

Datasheet

MYD88 purified MaxPab mouse polyclonal antibody (B01P)

Catalog Number: H00004615-B01P

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human MYD88 protein.

Immunogen: MYD88 (AAH13589, 1 a.a. ~ 296 a.a)
full-length human protein.

Sequence:

MAAGGPGAGSAAPVSSTSSLPLAALNMRVRRRLSLFL
NVRTQVAADWTALAEEMDFEYLEIRQLETQADPTGRL
LDAWQGRPGASVGRLLLELLTKLGRDDVLLELGPSIEE
DCQKYILKQQQEEAEKPLQVAADVSSVPRTAELAGITT
LDDPLGHMPERFDFICYCPSDIQFVQEMIRQLEQTN
YRLKLCVSDRDVLPGTCVWSIASSELIEKRCRRMVVVVS
DDYLSKECDFQTKFALSLSPGAHQKRLIPIKYKAMKK
EFPSILRSITVCDYTNPCTKSWFWTRLAKALSLP

Host: Mouse

Reactivity: Human, Rat

Applications: Det Ab, Flow Cyt, WB-Ce, WB-Ti, WB-Tr
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 4615

Gene Symbol: MYD88

Gene Alias: MYD88D