

Product Specification Sheet

VWR MICRO-WEIGHING DISHES, OVAL FORM

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

Cat. No.	Capacity (ml)	WxDxH (mm)	Pk
611-1357	0.04	6.5x3.5x2.5	250
611-1358	0.11	10x4x4.5	250
611-1360	0.12	12x4x2.7	250
611-1361	0.12	12x2x4.5	250
611-1359	1.12	25x7x7	50

Material: Aluminium

Other Characteristics:

- Ideal for microbalances, organic micro elementary and thermal analysis, and pharmaceutical use.
- These micro-weighing boats are compatible with most commercially available analysers and microbalances.
- Packaged in easy to use plastic containers for storage and dispensing.

Disclaimer: VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there is any doubt ask us.

Date of issue: **October 2013**