



Certificate of Conformity

VWR CRYOGENIC VIALS 2.0ml, SELF-STANDING, EXTERNAL THREAD POLYPROPYLENE

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

VWR Catalog No	Lot No
479-1208	225432088

Other Characteristics: **RNAse, DNAse, DNA and pyrogen free**

Minimum Specified Dose: **6.5 kGy**

Maximum Specified Dose: **12.5 kGy**

Quality Control:

Each shipment that goes to the sterilisation is verified by inserting a sample in a few cases of product to be shipped. These samples are then tested in our laboratory to verify their sterility.

All products that are then shipped to you have passed this sterilisation test. The Specific Dose of irradiation for this is a minimum of 6.5 and a maximum of 12.5kGy, for an average dose of 9.25kGy

This is to certify that this product is subject to proper manufacturing practices and inspection in accordance with the product specification using the sampling procedures MIL-STD-105E. It is inspected by qualified quality inspectors, and released after inspection and testing of determined criteria.

This document has been produced electronically and is valid without a signature

Date of issue: **December 2013**