

Product Specification Sheet

VWR CLEANROOM WIPES, NON WOVEN, SPEC-WIPE® 3

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

Cat. No.	Description	Sterility Level	LxW (mm)	Packed	Pk
115-0031	Spec-Wipe® 3	Non-sterile	229×229	300/bag	300
115-0032	Spec-Wipe® 3	Non-sterile	305×305	150/bag	150
115-0034	Spec-Wipe® 3, saturated with 70% IPA, 30% deionised water	Non-sterile	229×229	75/bag	75

Model	Spec-Wipe® 3
Absorbency - intrinsic (ml/g)	4.82
Absorbency - extrinsic (ml/m²)	335
Absorbency - rate (s)	1
Particle generation - 0.5-5 µm (×10 ⁶ /m²)	17.1
Extractables - using deionised water (g/m²)	0.056
Extractables - using isopropyl alcohol (g/m²)	0.015
Fibres (>100 µm) (×10 ³ /m²)	14.8
Ions - chloride (ppm)	7.8
Ions - sodium (ppm)	11.2
Basis weight (g/m²)	67

Product Specification Sheet

Other Characteristics:

- Suitable for use in Class 100, ISO 5 cleanrooms.
- Designed for cleaning laminar flow benches, cleanroom surfaces, quartzware and general process equipment.
- Appropriate for electronics, pharmaceutical, biotechnology, chemical and petroleum, environmental and education industries. Non textured.
- 45% polyester, 55% cellulose. Re-sealable pouches maintain product integrity, preserve cleanliness and saturation levels.
- Option of pre-saturated wipes, which reduce solvent use, chemical blending and waste.
- Hydroentangled without using chemicals, reduces ionic contamination

Disclaimer: VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there is any doubt ask us.

Date of issue: **December 2013**