



Material	24387.554
Material description	Glycerol
Grade	GPR RECTAPUR
CAS Number	56-81-5
Molecular formula	C ₃ H ₈ O ₃
Molecular mass	92.09

Characteristics	Specifications
Assay	Min. 98 %
Density (20/4)	1.257 --> 1.265
Heavy metals (as Pb)	Max. 20 ppm
Ignition residue (SO ₄)	Max. 0.05 %
Fe (Iron)	Max. 10 ppm

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Anja Vanhalle, Head of Laboratory - Haasrode
 VWR International bvba; Geldenaaksebaan 464; BE-3001
 Leuven; Belgium



For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.