



Material	84599.260
Material description	o-Cresolphthalein 20 g/l in ethanol 70%
Grade	VOLUSOL

Characteristics	Specifications
-----------------	----------------

Transition range pH 8.2 - 9.8

Passes test

## Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Anja Vanhalle, Head of Laboratory - Haasrode  
VWR International bvba; Geldenaaksebaan 464; BE-3001  
Leuven; Belgium

Transition range: pH 8,2 - 9,8

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.