



Material	95097.5000
Material description	Boric acid 4 % + indicator
Grade	VOLUSOL

Characteristics	Specifications
Assay (w/v)	3.9 --> 4.2 %
Appearance	Clear red liquid
Sensitivity	Passes test
Density (20/4)	1.000 --> 1.020

## Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Anja Vanhalle, Head of Laboratory - Haasrode  
VWR International bvba; Geldenaaksebaan 464; BE-3001  
Leuven; Belgium

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.