

FEP BOTTLE - 2L - 38-430



SPECIFICATIONS

UN-certification : 4G/X4.3/S/**/USA/+AA5031
4G/Y6.1/S/**/USA/+AA5031

Material: FEP

Neck Type : 38-430

Height : 246,00 mm

Length : 121,00 mm

Width : 121,00 mm

AVAILABLE ACCESSORIES :

Remark : Packaging materials and specifications shown may be modified without preceding notice
Version dd. 22 mai 2012

tel/tél +32(0)16 385 011
fax +32(0)16 385 385
info@be.vwr.com
www.vwr.com

Researchpark Haasrode 2020
Geldenaaksebaan 464 B-3001 Leuven
RPR/RPM Leuven
BTW/TVA BE 0403.593.343