



Certificate of Conformity

VWR CLEANING BRUSHES FOR TEST TUBES, WHITE BRISTLES WITH WIRE HANDLE

VWR International certifies that the products specified conform to all specifications and quality criteria mentioned in this document.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité mentionnés dans ce document.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle in diesem Dokument erwähnten Spezifikationen und Qualitätskriterien erfüllen.

Catalog No.	Material	Head Ø (mm)	Trimming length (mm)	Total length (mm)	Pk
149-0218	Natural bristle	13	75	200	5
149-0213	Natural bristle	19	81	230	5
149-0214	Natural bristle	35	100	250	5
149-0208	White nylon bristles	11	75	230	5
149-0219	White nylon bristles	13	78	200	5
149-0220	White nylon bristles	19	81	230	5
149-0221	White nylon bristles	35	100	250	5

Other Characteristics:

- Ideal for cleaning test tubes and a wide range of other laboratory utensils.
- With one or two tip ties to minimise risk of scratching and improve cleaning intensity
- Wire handle for secure grip
- With hanging loop

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **February 2015**