



Certificate of Conformity

VWR HUMIDIFIED SURFACE SAMPLING SPONGE ON HANDLE WITH NEUTRALIZING BUFFER IN STERILE SAMPLING BAG

VWR International certifies that the products specified conform to all specifications and quality criteria mentioned in this document.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité mentionnés dans ce document.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle in diesem Dokument erwähnten Spezifikationen und Qualitätskriterien erfüllen.

Cat. No.: 300-0232

Closure: 2 Round wire closure system

Closure Material: Steel wires covered with FDA approved pressure sensitive tape with a backing of polypropylene or PVC.

Dimensions and Volume (mm): 13.97 x 30.5 (cm)

Thickness of Material : 3 mil

Volume of Solution (Neutralizing buffer) : 10ml

Compliance: Comply to the applicable regulations of the Food and Drug Administration (21 CFR 177.1520 (c)3.1a, (c)3.2a and 178.2010) for use in contact with food.

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **August 2015**