



Certificate of Conformity

VWR PET SAMPLING BOTTLES, WIDE MOUTH WITH TAMPER EVIDENT SCREW CAP

VWR International certifies that the products specified conform to all specifications and quality criteria mentioned in this document.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité mentionnés dans ce document.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle in diesem Dokument erwähnten Spezifikationen und Qualitätskriterien erfüllen.

Cat. No.	Capacity (ml)	Dims HxWxL	Cap colour	Form	Pk	
331-0161	250	140X48X48	Blue	Square	200	
331-0155	500	172X62X62	Blue	Square	100	
331-0167	1000	243X73X73	Blue	Square	70	
331-0165	250	140X55	Orange	Round	200	
331-0166	500	168X70	Orange	Round	100	

Material: PET, Clear for Bottle, PP for caps

Properties:

- For basic sample collecting - Non sterile, lightweight, round or square.
- 34 mm neck for all bottles
- Tamper evident screw caps
- Leakproof
- **Do not autoclave**
- **Temperature resistance -10 to +40 Degrees C**

Other characteristics

- Clear PET ensures good visibility of the content level
- Ideal for water sampling and general sample collection
- Lightweight and durable
- Caps loose in bag in carton



Certificate of Conformity

Storage conditions :Room temperature. Do not expose to direct sunlight.

Autoclavable :No

Expiry :Not applicable

Product compliance:

Absence of heavy metals: Raw material is compliant with 94/62/EC directive

Food contact: Raw material is compliant to Regulation EC N°1935/2004

Standard NF en 13972: In rigid plastic containers. Nominal capacity definition, brim and total and measurement capabilities brim and total.

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **March 2015**