



Certificate of Conformity

VWR MEASURING BEAKERS WITH HANDLE

VWR International certifies that the products specified conform to all specifications and quality criteria mentioned in this document.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité mentionnés dans ce document.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle in diesem Dokument erwähnten Spezifikationen und Qualitätskriterien erfüllen.

Catalog No.	Capacity (ml)	D (mm)	H (mm)	Pk
213-0408	50	45	58	1
213-0409	100	71	72	1
213-0410	250	82	96	1
213-0411	500	97	131	1
213-0412	1000	127	160	1

Intended uses: Laboratory applications

Material: Porcelain C110

Shape: Round, with spout and handle

Compliant to: ISO material for the manufacture lies within the group C110 according to CSN EN 60672 –3
Product dimension tolerances are controlled by DIN 40680
Acid resistance acc. to ISO 1775

Properties:

Unglazed products can be used at temperatures up to 1350°C / 2462°F and glazed ones at temperatures up to 1000°C / 1.832°F

Highly resistant against acids and alkali

Density 2.4 g/cm³

Dending strength of unglazed products (70MPa)

Thermal conductivity (W.m⁻¹.K⁻¹) = 1.3



Certificate of Conformity

The compliance of demands is proved by regular tests and followed noticed tests from an independent laboratory:

Test according ČSN EN 60672-2
Test resistance to acids and bases

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **October 2015**