



Material	83626.320
Material description	Chloroform
Grade	HiPerSolv CHROMANORM for HPLC
CAS Number	67-66-3
Molecular formula	CHCl ₃
Molecular mass	119.38

Characteristics	Specifications
Assay (GC)	Min. 99.80 %
Water	Max. 0.0500 %
Non-volatile residue	Max. 0.001 %
Acidity	Max. 0.0005 meq/g
Alkalinity	Max. 0.0002 meq/g
Transmittance (250 nm)	Min. 50.0 %
Transmittance (260 nm)	Min. 84.0 %
Transmittance (280 nm)	Min. 95.0 %
Transmittance (300 nm)	Min. 98.0 %
Conforms to BDH 15283	Passes test

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

VWR International