



Material	87153.0500
Material description	Chloroform-[D1] 99.80 for NMR spectroscopy
CAS Number	865-49-6
Molecular formula	CDCl ₃
Molecular mass	120.38

Characteristics	Specifications
Assay (on anhydrous substance)	Min. 99.9 %
Isotopic enrichment (FT NMR 400 MHz)(D)	Min. 99.80 %
Water	Max. 0.01 %

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

VWR International