



Certificate of Conformity

VWR CULTURE TUBES WITH SCREW CAP

VWR International certifies that the products specified conform to all specifications and quality criteria mentioned in this document.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité mentionnés dans ce document.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle in diesem Dokument erwähnten Spezifikationen und Qualitätskriterien erfüllen.

VWR Catalog No.	ØxL(mm)	Thread	Pk
212-0656	12x100	GL 14	100
212-0657	16x100	GL 18	100
212-0658	16x160	GL 18	100
212-0659	18x180	GL 18	100

Material : Soda-lime glass, clear

Other Characteristics:

- Economical culture tubes.
- For use in applications without high thermal loads.
- Screw cap made from PP, with foam insert. Not acid resistant.

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **March 2016**