



## Certificate of Conformity

# VWR LDPE CARBOYS, NARROW MOUTH, ROUND, WITH SCREW CAPS, NON-STERILE, THREAD 83mm, THICKNESS 2.5 mm

VWR International certifies that the products specified conform to all specifications and quality criteria mentioned in this document.

*VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité mentionnés dans ce document.*

VWR International bestätigt hiermit, daß die angegebenen Produkte alle in diesem Dokument erwähnten Spezifikationen und Qualitätskriterien erfüllen.

Catalog No.	Capacity (L)	Neck I-Ø (mm)	Product Size (±2mm)	Height without Cap (±2mm)	Pk
216-1704	10	65	Ø250 X 385	379	1
216-1706	20	65	Ø287 X 535	520	1
216-1738	50	65	Ø378 X 680	663	1

### Material

**Carboy** : LDPE, 100 % Virgin, FDA 21 CFR Compliance, USP Class VI passed.

**Closure** : Polypropylene, 100 % Virgin, FDA 21 CFR Compliance, USP Class VI passed With Food grade colour.

**BSE/TSE Statement:** These product is manufactured without the use of Raw materials of animal origin.

**Latex Free Statement:** No material containing latex or natural rubber are used in the manufacturing, handling and packaging process for this product.

### Other Characteristics:

- Preferred for high purity applications as it has high resistance to chemicals (eg. widely used in ophthalmic product range, like eye droppers, nasal spray, etc.)
- Preferred resin for cold storage (up to -80 °C), working temperature +80 °C to -80 °C
- Resistance to physical impacts ( products will not break even if it falls down by accidentally, where as PP carboys may lead to crack at negative temperature)
- LDPE product will bounce back to its shape as its high flexible material
- Widely used in pharma industries

**Disclaimer:** VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **December 2015**