



K-400 P

Model	K-400 P
Optical System	Common main objective (CMO) Infinity optical system
Observation Tube	Binocular head
Interpupillary distance	54-76mm
Diopter adjustment	On both tubes , +/- 5 diopter
Inclination	45° inclined
Eyepieces	WF10X/23mm
Objectives system	4-position step magnification changer
Standard objective	1X Achromatic
Objectives magnification	6X, 12X, 25X and 50X
Working distance	89mm
Stand	Pole type stand
Base	329x277mm
Column/Arm	Ø32mm, 295mm high
Optical body with head holder	For Ø32mm column
Focus mechanism	Coarse focusing system with tension adjustment
Focusing stroke	49mm
Accessories included	Black and white stage plates and dust cover
Length x Width x Height	329x277x399mm
Weight	5,8Kg