



Certificate of Conformity

VWR DOUBLE POUCH TWIST SEAL BAG, STERILE, SAFETY TABS, PRINTED, CLOSURE WITH 1 ROUND WIRE AND 1 FLAT WIRE

VWR International certifies that the products specified conform to all specifications and quality criteria mentioned in this document.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité mentionnés dans ce document.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle in diesem Dokument erwähnten Spezifikationen und Qualitätskriterien erfüllen.

Cat. No.	Capacity (ml)	Thickness (µm)	Dimensions (mm)	Pk
129-9012	650	76	140x229	1000
129-9013	650	101	140x229	1000
129-9014	1080	76	140x305	1000
129-9015	1650	76	178x305	1000
129-9016	1650	101	178x305	1000
129-9017	2250	101	254x305	1000
129-9018	4000	101	254x305	1000

Bag material : Polyethylene, transparent

Sterility :

Product is gamma irradiated at a min level of 5 Kgy. Products is packaged in Case of 1000 bags (40 packs of 25 bags) and each pack is double pouched, divided into 2 boxes (20 packs per box) .Five years sterility guaranteed.

Other characteristics:

- Offer an economical and efficient way to collect , contain and carry samples.
- Maximum bag mouth opening to facilitate insertion of samples
- Unique fold over system providing a temporary leak-proof seal
- Bags accommodate solid, semi solid and liquid
- Sterile material
- R-Nase, D-Nase and Pyrogene free

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **July 2016**