



Certificate of Conformity

VWR MAGNETIC STIRRING BARS, CROSS SHAPED

VWR International certifies that the products comply in every respect with the relevant particular specifications and/or drawings and/or relevant standards and regulations in force, mentioned in this document.

VWR International certifie que les produits sont conformes à tous égards aux spécifications particulières appropriées et / ou dessins, normes et règlements en vigueur, mentionnées dans le présent document

VWR International bestätigt, dass das Produkte allen zutreffenden Spezifikations- und/oder Zeichnungs- und/oder Standardanforderungen sowie den im Dokument aufgeführten gültigen Vorschriften entspricht.

VWR Cat.No.	Height (mm)	Ø (mm)	Pk
442-0394	5	10	10
442-0395	8	20	10
442-0396	9	25	10
442-0397	10	30	10
442-0398	11	38	10
442-0399	15	50	10
442-0400	20	60	10

Material : PTFE*, conforming to FDA and USP Class VI requirements

Sterilisation: Chemical and thermal NOT gamma radiation

Other characteristics :

- Smooth surface.
- Stable rotating position for optimal stirring effect.
- Very good chemical resistance.
- Temperature resistant, suitable for high and low temperature applications (–200 to +280 °C).

* PTFE used in the manufacture of this product is VIRGIN PTFE GRADE TF1620 manufactured by 3M® DYNELON®

Certificate of Conformity

3M® DYNEON® TF1620 – Compliance Status

FDA – Directive FDA 21 CFR 177.1550

USP Class VI – USP 27, NF 22, 2004 for Plastics at 70°C.
RoHS 2 – Directive 2011/65/EU & Dangerous Substances – 76/769/EC
EU REACH Regulations – SVHC Current Candidate List
Bisphenol A – Absent
Latex – Absent
ADCF – Animal Derived Components Free (BSE/TSE). GMO Free
Endotoxins – Free
POP, PFOS, PFOA, APFO, DEHIP, VOC Free
ODC – Class 1
Dodd-Frank Act – Gold, Coltan, Cassiterite, Wolframite and their derivatives originating from the Democratic Republic of Congo - Absent

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **September 2016**