	Document Number	Revision	Page	Date
	CSP001	A0	1 of 1	12-04-2016
PRODUCT DATA SHEET		Cell Strainer Pestle		


1: Material:

Pestle: Polypropylene PP Color: Purple(259C)

2: Related Catalog Numbers:

VWR EU EAN. No.	Description	Sterile	Packout Description
734-2819	Pestle for Cell Strainer, Purple, Sterile, Individually Packaged	Yes	Individually Packaged, 100/case
734-2820	Pestle for Cell Strainer, Purple, Sterile, Bulk	Yes	25/bag, 50/case

3: Dimensions=mm



4: Features

- Convex Pestle head includes molded textured surface to better manipulate material in the cell strainer
- Made from high grade gamma resistance polypropylene
- Autoclavable at 121 °C
- Individually packed in peel-to-open paper/plastic wrapper and bulk pack with ziplock bag
- 5.5 in length for convenient hand manipulation
- Certified DNase/RNase Free and Non-pyrogenic
- Shelf life: 3 years after month of production

5: STERILITY

- Sterilized by gamma irradiation SAL10⁻⁶ (ISO11137)

6: PRODUCT ISO 9000 CERTIFICATION

- Manufactured in a class 100,000 room environment.
- Manufactured under EN ISO13485 :2003 and ISO9001:2008 quality management system.

Drawn	Approved
Date 12 April 2016	Date 12 April 2016
Name Alex Fang	Name James