

Certificate of Conformity

BOUFFANT CAP LATEX FREE 61 CM 113-7919

CHARACTERISTICS:

Bouffant clip caps, are made of 11.5 grams/m² by weight, non-woven polypropylene material, commonly used in a clean environment. Bouffant clip caps are manufactured in a stringent clean environment. Color white.

PHYSICAL PROPERTIES

Material	Non-woven Polypropylene 11.5grams	
Elastic	latex free	
Size	61 cm + /- 1 (Length), 23 cm +/- 1 (width)	
Stretch Length	61 cm +/- 1	
Colors	White	
Tensile Strength	MD > 16N/5cm;	CD > 10N/5cm
Ultimate Elongation	>30%	

PACKING:

100pcs/bag, 10 bags/case

QUALITY ASSURANCE CRITERIA:

Bouffant caps are guaranteed to meet 4.0 AQL level II sampling plan per MIL-std-150E for the following defect criteria:

- | | |
|--------------------------|---|
| 1. Broken elastic band | 4. Holes/tears in PP material caused by machine |
| 2. Irregular PP material | 5. Physical dimensions non-compliance |
| 3. Dirty PP material | |

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **November 2015**