



Certificate of Conformity

VWR BEAKERS, LOW FORM, PTFE, WITH SPOUT

VWR International certifies that the products comply in every respect with the relevant particular specifications and/or drawings and/or relevant standards and regulations in force, mentioned in this document.

VWR International certifie que les produits sont conformes à tous égards aux spécifications particulières appropriées et / ou dessins, normes et règlements en vigueur, mentionnées dans le présent document

VWR International bestätigt, dass das Produkte allen zutreffenden Spezifikations- und/oder Zeichnungs- und/oder Standardanforderungen sowie den im Dokument aufgeführten gültigen Vorschriften entspricht.

Cat. No.	Capacity (ml)	Height (mm)	Ø ext. (mm)	Pk
213-0233	5	26	20	1
213-0234	10	33	24	1
213-0235	25	47	32	1
213-0236	50	60	43	1
213-0237	100	68	54	1
213-0238	150	69	59	1
213-0239	250	97	66	1
213-0240	500	125	80	1
213-0241	600	125	90	1
213-0242	1000	155	100	1
213-0243	2000	210	120	1
213-0244	5000	280	170	1
213-02700	1	18	13	2
213-0274	400	106	80	1

Material: PTFE conforming to FDA and USP Class VI requirements

Minimum / maximum use Temp.: From -200 up to 300°C

Other characteristics:

- Excellent chemical resistance.
- Temperature resistant from -200 to +280 °C.
- Flat bottom for good heat transfer, for example when used on regulated hotplates.



Certificate of Conformity

PTFE used in the manufacture these products are pure, virgin material which conforms to FDA and the following specified substance content. It does not contain, or use in its manufacture any Class 1 ODC according to the Clean Air Act 1990, PFOA, POP, DEHIP and BPA, no animal material is used in the manufacture or is present in this material.

	Description	Content (%)
Group A	Cadmium and Cadmium compounds	Not detected
	Hexavalent Chromium compounds	Not detected
	Lead and Lead compounds	Not detected
	Mercury and Mercury compounds	Not detected
	Tributyl Tin Oxide (TBTO)	Not detected
	Tributyl Tins & Triphenyl Tins	Not detected
Group B	Polybrominated Biphenyls (PBBs)	Not detected
	Polybrominated Diphenyl ethers (PBDE)	Not detected
	Polychlorinated Biphenyls (PCBs)	Not detected
	Polychloronaphthalenes	Not detected
	Short Chain Chlorinated Paraffins	Not detected
Group C	Asbestos	Not detected
	Azo Colorants	Not detected
	Ozone Depleting Substances	Not detected
	Radioactive Substances	Not detected

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **January 2018**