



Certificate of Conformity

VWR CONICAL SAMPLE CONTAINERS WITH SNAP-ON LIDS

VWR International certifies that the products specified conform to all specifications and quality criteria mentioned in this document.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité mentionnés dans ce document.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle in diesem Dokument erwähnten Spezifikationen und Qualitätskriterien erfüllen.

Catalog No.	Capacity (ml)	Height (mm)	Type	Ø top/ Ø bottom (mm)	Packed	Pk
216-2648	200	88	Separate snap-on lid	62/54	Carton of 660	660
216-2649	200	88	Assembled snap-on lid, sterile	62/54	Carton of 220	220
216-2651	400	100	Separate snap-on lid	85/72	Carton of 460	460
216-2652	400	100	Assembled snap-on lid, sterile	85/72	Carton of 185	185
216-2639	1000	130	Separate snap-on lid	105/96	Carton of 250	250
216-2640	1000	130	Assembled snap-on lid, sterile	105/96	Carton of 80	80

Container Material: **Natural polypropylene**

Cap Material: **LDPE**

Physical properties: **Temperature resistance main body -20 °C / + 120 °C**
Temperature resistance secondary body -80 °C / + 75 °C
Flexible to rigid, good mechanical resistance, high chemical resistance
Main body is autoclavable

Compliant to: **Raw material is compliant with the EC/94/62 directive on absence of heavy metals**
Raw material is compliant with regulation EC N°1935/2004 on food contact
NF EN 13972 nominal capacity definition, brim and total measurement capabilities brim and total

If sterile, ISO 11137 Sterilisation of health care products

Shelf Life: **3 years**

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **February 2018**