



Material	443815S
Material description	D(+)-Sucrose
Grade	Electran Molecular biology grade
CAS Number	57-50-1
Molecular formula	C <sub>12</sub> H <sub>22</sub> O <sub>11</sub>
Molecular mass	342.30

Characteristics	Specifications
DNases (exo- and endonucleases)	Not detected
Identification (IR)	Conforms
Proteases	Not detected
RNases:	Not detected
Specific optical rotation (26 %; water)	66.3 --> 67.0 °
Fe (Iron)	Max. 0.1 ppm
Pb (Lead)	Max. 0.5 ppm

## Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

VWR International