

Certificate of Conformity

VWR Catalogue Number	See List
Description	Measuring Beakers with Handle, Embossed scale, Blue Graduations and Spout
Country of Origin	Manufactured in Germany
Date of Issue (yyyy-mm-dd)	2021-08-03

Catalogue Number	Capacity [mL]	Height [mm]	External Diameter [mm]	Graduation Interval [mL]
213-3401	250	120	70	5
213-3402	500	133	91	10
213-3403	1000	182	116	20
213-3404	2000	215	150	50

Quality System Compliance

Products are manufactured under the **ISO 9001:2015** certified Laboratory. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material

Polypropylene (PP); highly translucent.

The material used to produce the products are listed in List of Authorised Substances of Regulation EU/10/2011.

The material grades conform to EU Directive 1935/2004 on materials and articles intended to come into contact with food.

Temperature Resistance

+121 °C (for 20 minutes) prolonged use at +100 °C.

Sterilization	Not applicable.
ATP Assay	Not applicable.
DNase & RNase Free	Not applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Not applicable.
Latex Statement	Not applicable.
BPA Statement	Not applicable.
DEHP Statement	Not applicable.
RoHS	Not applicable.
REACH Statement	Not applicable.
Storage Conditions	Depends on usage.
Shelf Life	Depends on usage.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature