



Certificate of Conformity

VWR STERILE DISPOSABLE SPATULAS

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

VWR Catalog No:	Lot No:	Colour	Length
231-0369	122011	green	210 mm

Material: **100% Food grade Virgin Polypropylene**
Raw materials free from any animal origin

Compliant to: **FDA standards for plastics coming in contact with food**

Quality control for the product is performed upon receipt of each batch. Random samples of each batch are tested for adherence to size, weight and other manufacturing specifications; ability to withstand boiling, autoclaving and exposure to liquid nitrogen without damage to the product are also tested.

This document has been produced electronically and is valid without a signature

Date of issue: **April 2012**