



Product Specification Sheet

VWR CULTURE AND SAMPLE TUBES, STERILE, WITHOUT CAPS

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

VWR Cat.No.: **211-0133**

Product Dimensions: **Length of Tube: 95.65mm**
Outer Diameter: 18.00mm
Inner Diameter: 15.65mm

Packaging: **125 tubes per bag**
8 bags per case (1000 tubes total)
Can be centrifuged to 1500 RCF

Sterilization: **Sterility assurance level (SAL) between 10^{-5} and 10^{-6} for all our sterile products.**

Pyrogens: **Certified testing is performed lot by lot. These tests include Nuclease(Rnase Dnase) and Endotoxin (pyrogen). Certificates are on file at the production facility.**

Performance Testing: **Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures.**

Disclaimer: *VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there is any doubt ask us.*

Date of issue: **June 2012**



VWR International bvba/sprl, Researchpark Haasrode 2020, Geldenaaksebaan 464, 3001 Leuven, Belgium
<http://www.vwr.com>

Version: 2.0