

VWR®

Product Specifications Sheet

VWR WOMEN'S CHEM COAT WITH STAND-UP COLLAR

Coat women's style with stand-up collar made from 99% Polyester fibres with 1% Hydro-Tec finish (carbon treated), 2 side-pockets + 1 breast-pocket on the left (all covered by a flap), long sleeves with elastic cuffs, concealed press stud band with 5 V2A-steel press studs (covered by a flap), fitted, with back belt

Intended use

General laboratory uses
Protection against liquid chemical splashes
Chemical manufacturing areas

Specifications

Material	99% Polyester with 1% Hydro-Tec finish EN 13034:2005 (Type 6), EN 1149-5
Weight (approx.)	275g/m ²
Color	White
Fibre characteristic	OEKO-TEX Standard 100
Washing instruction	 suitable for industrial laundering



VWR®

Product Specifications Sheet

VWR WOMEN'S CHEM COAT WITH STAND-UP COLLAR

Sizes				Measurements in cm				Catalogue number
D	IT	UK	FR	½ Chest size (1)	Crown width (2)	Total (3)	Arm length (incl. cuffs) (4)	
34	40	8	36	47,5	37,9	95,5	57	113-1825
36	42	10	38	49,5	39,3	97	58	113-1826
38	44	12	40	51,5	40,7	98,5	59	113-1827
40	46	14	42	53,5	42,1	100	60	113-1828
42	48	16	44	55,5	43,5	101,5	61	113-1829
44	50	18	46	57,5	44,9	103	62	113-1830
46	52	20	48	59,5	46,3	104,5	63	113-1831
48	54	22	50	61,5	47,7	106	64	113-1832
50	56	24	52	63,5	49,1	107,5	65	113-1833
52	58	26	54	65,5	50,5	109	66	113-1834



EN 13034:2005+A1:2009

Typ PB [6]



EN 1149-5:2008



Product compliance:

Council Directive 89/686/EEC relating to Personal protective equipment, CAT III against high risks, EN 13034 Type 6: Partially splash-proof chemical protective clothing
Certification on Personal Protective Equipment + EN13034 available.

Disclaimer: VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there is any doubt ask us.