

## Bolton broth selective supplement

**Code 84744.0001**

### Intended use

Freeze-dried selective supplement for the preparation of Bolton broth used for the enrichment of *Campylobacter* spp. in food and animal feeding stuffs (ISO10 272).

### Vial contents

Cefoperazone.....10.0 mg  
Vancomycin..... 10.0 mg  
Trimethprim lactate .....10.0 mg  
Amphotericin B.....5.0 mg

### Instructions for preparation

Aseptically reconstitute the contents of one vial with 5 ml of ethanol/sterile purified water (1:1). Mix well and aseptically add to 500 ml of Bolton broth base Art.: 84696.0500, autoclaved and cooled to approximately 45-50 °C. Add also 25 ml of lysed horse blood. Mix well and pour sterile tubes or flasks under aseptic conditions.

### Instructions for use

For laboratory use only.

Consult the technical sheet of Bolton broth base Art.: 84696.0500

Consult the material safety data sheet before the use.

### Physical characteristics

|   |   |
|---|---|
| Appearance of freeze-dried supplement       | Short, pale yellow pastille               |
| Appearance of solution after reconstitution | Yellow, opalescent with white precipitate |
| pH of complete medium (25°C)                | 7.4 ± 0.2                                 |

### Storage conditions

Store at 2-8°C. When stored as directed the supplement remains stable until the expiry date shown on the label. Do not use beyond stated expiry date.

### Ordering information

84744.0001    Bolton broth selective supplement    10 vials, each for 500 ml of complete medium