

## Datasheet

### Sfpi1 polyclonal antibody

**Catalog Number:** PAB5935

**Regulation Status:** For research use only (RUO)

**Product Description:** Rabbit polyclonal antibody raised against synthetic peptide of Sfpi1.

**Immunogen:** A synthetic peptide corresponding to N-terminus of mouse Sfpi1.

**Host:** Rabbit

**Reactivity:** Mouse

**Applications:** ELISA, WB

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Recommend Usage:** The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS, pH 7.2 (50% glycerol, 0.01% sodium azide)

**Storage Instruction:** Store at -20°C.  
Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 20375

**Gene Symbol:** Sfpi1

**Gene Alias:** Dis-1, Dis1, PU.1, Sfpi-1, Spi-1, Tcfpu1, Tfpu.1