

Datasheet

AGR2 polyclonal antibody (A01)

Catalog Number: H00010551-A01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a partial recombinant AGR2.

Immunogen: AGR2 (NP_006399, 21 a.a. ~ 110 a.a) partial recombinant protein with GST tag.

Sequence:

RDTTVKPGAKKDTKDSRPKLPQTLSRGWGDQLIWTQ
TYEEALYKSKTSNKPLMIHHLDECPHSQALKKVFAEN
KEIQKLAEQFVLLNLV

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Re

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: 50 % glycerol

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 10551

Gene Symbol: AGR2

Gene Alias: AG2, GOB-4, HAG-2, XAG-2