

## Datasheet

### HTATIP2 monoclonal antibody (M01), clone 4G8

**Catalog Number:** H00010553-M01

**Regulation Status:** For research use only (RUO)

**Product Description:** Mouse monoclonal antibody  
raised against a full length recombinant HTATIP2.

**Clone Name:** 4G8

**Immunogen:** HTATIP2 (AAH02439, 1 a.a. ~ 242 a.a.)  
full-length recombinant protein with GST tag. MW of the  
GST tag alone is 26 KDa.

**Sequence:**

MAETEALSKLREDFRMQNKSVFILGASGETGRVLLKEI  
LEQGLFSKVTLIGRRKLTDFDEEAYKNVNQEVVDFEKLD  
DYASAFQGHVDVGFCCLTTRGKAGAEFVRVDRDYV  
LKSAELAKAGGCKHFNLLSSKGADKSSNFLYLQVKGE  
VEAKVEELKFDRYSVFRPGVLLCDRQESRPGEWLVRK  
FFGSLPDSWARGHSVPVTVVRAMLNNVVRPRDKQM  
ELLENKAIHDLGKAHGSLKP

**Host:** Mouse

**Reactivity:** Human

**Applications:** ELISA

(See our web site product page for detailed applications  
information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product  
page for detailed protocols

**Isotype:** IgG2a lambda

**Storage Buffer:** In 1x PBS, pH 7.4

**Storage Instruction:** Store at -20°C or lower. Aliquot to  
avoid repeated freezing and thawing.

**Entrez GeneID:** 10553

**Gene Symbol:** HTATIP2

**Gene Alias:** CC3, FLJ26963, SDR44U1, TIP30