

Datasheet

ZMYM3 polyclonal antibody

ZNF198L2, ZNF261

Catalog Number: PAB19387**Regulation Status:** For research use only (RUO)**Product Description:** Rabbit polyclonal antibody raised against synthetic peptide of ZMYM3.**Immunogen:** A synthetic peptide corresponding to 20 amino acids near C-terminus of human ZMYM3.**Host:** Rabbit**Reactivity:** Human, Mouse**Applications:** ELISA, WB-Ti

(See our web site product page for detailed applications information)

Protocols: See our web site at<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols**Specificity:** At least three isoforms of ZMYM3 are known to exist; this antibody will detect all three. ZMYM3 antibody is predicted to not cross-react with other ZMYM protein family members.**Form:** Liquid**Purification:** Peptide affinity purification**Concentration:** 1 mg/mL**Recommend Usage:** Western Blot (1 ug/mL)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS (0.02% sodium azide)**Storage Instruction:** Store at 4°C for three months. For long term storage store at -20°C.
Aliquot to avoid repeated freezing and thawing.**Entrez GeneID:** 9203**Gene Symbol:** ZMYM3**Gene Alias:** DXS6673E, KIAA0385, MYM, XFIM,