



Material	05547535
Material description	Isomount 2000
Grade	Q path® mounting medium

Characteristics	Specifications
Appearance	Colourless liquid
Viscosity (20°C) (day 0)	300 --> 400 mPa.s
Viscosity (20°C) (day 3)	400 --> 500 mPa.s

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Sewing
Quality Control - VWR International

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.