



Certificate of Conformity

VWR BENCH PROTECTORS

VWR International certifies that the products specified conform to all specifications and quality criteria mentioned in this document.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité mentionnés dans ce document.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle in diesem Dokument erwähnten Spezifikationen und Qualitätskriterien erfüllen.

Catalog No.	L×W (mm)	Packed	Pk
115-9221	50000×460	1 roll of 50 m	1
115-9220	570×460	50 single sheets	50

Colour : White

Absorbance rate: 500 ml/m²

Grammage, gsm : 165 - 185

Thickness, µm@53Kpa : < 350

Other Characteristics:

- High quality, smooth, absorbent paper that quickly absorbs liquid spills with a thick, laminated polyethylene layer that prevents seeping onto the work surface.
- Retains leaked reagents in radio-chemical laboratories and prevents radioactive contamination
- Enables recovery of leaked, expensive materials
- Breakage protection for hard surfaces
- Absorption layer for water or solvent moisture chambers
- Ideal, impermeable protective material for lining animal cages

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **March 2015**