



Certificate of Conformity

VWR WIPES FOR GENERAL APPLICATION

VWR International certifies that the products specified conform to all specifications and quality criteria mentioned in this document.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité mentionnés dans ce document.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle in diesem Dokument erwähnten Spezifikationen und Qualitätskriterien erfüllen.

Catalog No.	Description	L×W (mm)	Pk
115-0203	Tissue wipes, 1 ply	424×355	140
115-0202	Tissue wipes, 1 ply	213×114	280
115-0204	Tissue wipes, 3 ply	424×368	90

Material: 100% Virgin Fiber Paper

Colour: White

Other Characteristics:

- Absorbent, soft, thick and low-lint fleece wipes.
- For light use in the laboratory, electronics, computer industry, pharmaceuticals, biotechnology, optics, etc.
- High absorbency, in three formats, available in single or triple ply.
- Practical dispenser box with anti-static, plastic film window for removal.
- Table and wall-mounted dispensers for one handed removal of the wipes are available separately.

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **January 2016**