

## Datasheet

### Slc12a1 polyclonal antibody (ATTO 680)

**Catalog Number:** PAB29043

**Regulation Status:** For research use only (RUO)

**Product Description:** Rabbit polyclonal antibody raised against synthetic peptide of Slc12a1.

**Immunogen:** A synthetic peptide corresponding to amino acids 33-55 at N-terminus of rat Slc12a1.

**Host:** Rabbit

**Reactivity:** Human, Mouse, Rat

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Specificity:** Detects ~160 kDa.

**Form:** Liquid

**Conjugation:** ATTO 680

**Purification:** Affinity purification

**Recommend Usage:** Western blot (1:1000)  
The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS, pH 7.4 (50% glycerol, 0.09% sodium azide)

**Storage Instruction:** Store at 4°C.  
Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 25065

**Gene Symbol:** Slc12a1

**Gene Alias:** MGC93479, Nkcc2