

Datasheet

Slc12a1 polyclonal antibody (ATTO 700)

Catalog Number: PAB29044

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of Slc12a1.

Immunogen: A synthetic peptide corresponding to amino acids 33-55 at N-terminus of rat Slc12a1.

Host: Rabbit

Reactivity: Human, Mouse, Rat

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: Detects ~160 kDa.

Form: Liquid

Conjugation: ATTO 700

Purification: Affinity purification

Recommend Usage: Western blot (1:1000)
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (50% glycerol, 0.09% sodium azide)

Storage Instruction: Store at 4°C.
Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 25065

Gene Symbol: Slc12a1

Gene Alias: MGC93479, Nkcc2