

## Datasheet

### c-Myc polyclonal antibody

**Catalog Number:** PAB29060

**Regulation Status:** For research use only (RUO)

**Product Description:** Chicken polyclonal antibody raised against recombinant c-Myc peptide.

**Immunogen:** A KLH-conjugated synthetic peptide corresponding to c-Myc tag .

**Sequence:** EQKLISEEDL

**Host:** Chicken

**Applications:** ELISA, ICC, IHC, IP, WB  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at  
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Purification:** Antigen affinity purification

**Isotype:** IgY

**Recommend Usage:** ELISA

Immunocytochemistry(1:250)

Immunohistochemistry

Immunoprecipitation(1:250)

Western Blot(1:500)

The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS, pH 7.2 (0.02% sodium azide)

**Storage Instruction:** Store at 4°C and avoid from light.

For long term storage, dilute the antibody with sterile glycerol and store it at -20°C.

Aliquot to avoid repeated freezing and thawing.