



Certificate of Conformity

VWR CRYOGENIC VIALS WITH EXTERNAL THREAD AND LIP SEAL

VWR International certifies that the products specified conform to all specifications and quality criteria mentioned in this document.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité mentionnés dans ce document.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle in diesem Dokument erwähnten Spezifikationen und Qualitätskriterien erfüllen.

Catalog No.	Capacity (ml)	Base	Ø ext.×H (mm)	Pk
479-1219	1.2	Freestanding	12.5×41	1.000
479-1235	2.0	Round bottom	12.5×46	1.000
479-1220	2.0	Freestanding	12.5×48	1.000
479-1775	3.0	Freestanding	12.5×71	1.000
479-1776	4.0	Freestanding	12.5×77	1.000
479-1777	5.0	Round bottom	12.5×91	1.000
479-1236	5.0	Freestanding	12.5×91	1.000

Material: **PP vial and closure, sterile (SAL 10⁻⁶)**

Temperature range: **-196°C to +121°C Autoclavable at 121°C, for up to 30 minutes.**

Product compliance:

Resin conforms to FDA 21 CFR 177.1520

Resin conforms to USP Class VI

Plastics and colorants are in conformity with CONEG standards

Other Characteristics:

Certified RNase, DNase, Pyrogen and DNA free

Gamma radiation sterilized.

Printed graduations and writing area

Clear , no presence of contamination in plastic.

Closure verification

Leak proof testing in vacuum of 71.3cm/Hg.

Gas phase of Liquid nitrogen resistance testing.

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **June 2016**